

Premisut Surgical Latex Gloves

Description

Premisut Surgical Latex sterile gloves are essential protective gear for healthcare professionals during surgical procedures. Available in both pre-powdered and powder-free variants, these gloves are meticulously crafted from latex, ensuring flexibility, durability, and a snug fit. Pre-powdered gloves facilitate easy donning, while powder-free options eliminate the risk of powder-related complications. Individually packaged to maintain sterility, they provide a reliable barrier against contaminants, safeguarding both patient and practitioner from potential infections.

Premisut Surgical
Latex Gloves
(Pre Powdered)

Premisut Surgical
Latex Gloves
(Powder free)



Premisut SURGICAL LATEX GLOVES

Brand Name	PREMISUT SURGICAL GLOVES POWDERED						
Conformity	Surgical Rubber Gloves meets all the requirements of International std. and National std. ASTM D 3577, EN 455-1,2,3 & 4, ISO 10282 & IS 13422 standards.						
CE Certificate	CE 2265 & ISI Mark						
Classification	Premisut Surgical Gloves Powdered are classified as Class IIa device.						
Colour & Visual appearance	Pale Yellow Anatomical shape, Micro Roughened surface and High strength						
Physical dimension's	Size	Length (mm)	Width(mm)	Thickness (mm)		Bead thick	
	6.0	280 min.	77 +/- 5	Cuff	0.10 min.	1.50 – 2.00 mm	
	6.5		83 +/- 5	Palm	0.10 min.		
	7.0		89 +/- 5				
	7.5		95 +/- 5				
	8.0		102 +/- 5	Finger	0.10 min.		
	8.5		108 +/- 5				
	9.0		114 +/- 5				
Physical properties	Tensile strength MPa		Stress @ 500% MPa	Ultimate Elongation (%)		Force @Break (N)	
Before Ageing	24.0 min.		5.5 max.	750 min.		9.0 min.	
After Ageing	18.0 min.		--	560 min.		9.0 min.	
Other properties	Powder content		Protein content	Holes Level		Moisture	
	15.0 mg/dm2 max.		200 µg/dm2 max.	AQL-0.65 (EN 455-1) AQL- 1.5 (ASTM D 3577)		1.2% max.	
Glove average weight size 7.5	18.0 g/pair min.						
Sterilization Method	Ethylene Oxide						
Storage	Recommended storage condition: 0 - 40° C, keep away from moisture and direct heat						
Shelf life	5 years						
Packing	Wrapper packed one pair of gloves primarily packed in sterile pack then it is sealed, 50 unit of sterile packs are packed in a printed dispensary box. 08/10/12 such boxes are packed in a master carton – 400/500/600 Pairs/Shipper box.						
	Wrapper pack		Sterile pack		Dispensary pack		Shipment pack
	MG poster paper		Tear pack/ Peel pack		3 – Ply		5 – Ply
POWDER	SIZES						
Sterile Disposable Latex Surgical Gloves	6.0	6.5	7.0	7.5	8.0	8.5	9.0
Ref	AMSGP6.0	AMSGP6.5	AMSGP7.0	AMSGP7.5	AMSGP8.0	AMSGP8.5	AMSGP9.0

Premisut SURGICAL LATEX GLOVES

Brand Name	PREMISUT SURGICAL GLOVES POWDER FREE						
Conformity	Surgical Rubber Gloves meets all the requirements of International std. and National std. ASTM D 3577, EN 455-1,2,3 & 4 and ISO 10282 standard.						
CE Certificate	CE 2265						
Classification	Premisut Surgical Gloves Powder Free are classified as Class IIa device.						
Colour & Visual appearance	Pale Yellow Anatomical shape, Micro Roughened surface and High strength						
Physical dimension's	Size	Length (mm)	Width(mm)	Thickness (mm)		Bead thick	
	6.0	280 min.	77 +/- 5	Cuff	0.10 min.	1.50 – 2.00 mm	
	6.5		83 +/- 5	Palm	0.10 min.		
	7.0		89 +/- 5				
	7.5		95 +/- 5	Finger	0.10 min.		
	8.0		102 +/- 5				
	8.5		108 +/- 5				
	9.0		114 +/- 5				
Physical properties	Tensile strength MPa		Stress @ 500% MPa	Ultimate Elongation (%)		Force @Break (N)	
Before Ageing	24.0 min.		5.5 max.	750 min.		9.0 min.	
After Ageing	18.0 min.		--	560 min.		9.0 min.	
Other properties	Powder content		Protein content	Holes Level		Moisture	
	2.0 mg/glove max.		50 µg/dm2 max.	AQL-0.65 (EN 455-1) AQL- 1.5 (ASTM D 3577)		1.2% max.	
Glove average weight size 7.5	18.0 g/pair min.						
Sterilization Method	Ethylene Oxide						
Storage	Recommended storage condition: less than 40° C, keep away from moisture and direct heat						
Shelf life	5 years						
Packing	Wrapper packed one pair of gloves primarily packed in sterile pack then it is sealed, 50 unit of sterile packs are packed in a printed dispensary box. 8/12 such boxes are packed in a master carton – 400/600 Pairs / Shipper box.						
	Wrapper pack		Sterile pack		Dispensary pack		Shipment pack
	MG poster paper		Tear pack/ Peel pack		3 – Ply		5 – Ply
POWDER FREE				SIZES			
Sterile Disposable Latex Surgical Gloves	6.0	6.5	7.0	7.5	8.0	8.5	9.0
Ref	AMSGPF6.0	AMSGPF6.5	AMSGP7.0	AMSGPF7.5	AMSGPF8.0	AMSGPF8.5	AMSGPF9.0