

Premisut Bone Wax

Premisut

Description

Premisut BONE WAX is a sterile bone wax, mixture of White Bees wax, Hard Paraffin and Isopropyl Palmitate, a wax-softening agent. Premisut BONE WAX is opaque and has a waxy odor.

Indications

Premisut Bone Wax may be used for the control of bleeding from bone surfaces.

Performance

Premisut BONE WAX achieves local haemostasis of bone by acting as a mechanical (tamponade) barrier. It does not act bio chemically and is minimally resorbable.

Sterilization
Ethylene Oxide

Shelf Life
5 years

Packaging
Box of 12pcs

REF	SIZE	PCS/BOX
PBW-810	2.5gr	12

