


Description

Premisut ORC is a natural, plant based Absorbable Haemostat, made up of Oxidized Regenerated Cellulose designed to stop capillary, venous and small arteriolar bleeding, when ligation and other conventional procedures are either ineffective or impractical.

Features & Benefits:

- Rapid Actions: haemostasis within 3-4 min
- Absorbable: completely
- Biocompatible: Minimum tissue reaction
- Antibacterial: Low pH leads to the bacterial effect
- Adaptable: Can be cut into any desire dimensions without fraying and applicable to various tissue thickness and surgery sites
- Easy to use: Highly flexible and drapable; can be sutured
- Ready to use: Applied, laid on, held against or wrapped around the bleeding surface in dry form
- Safer: No animal origin or collagen substances
- Compliance: USP & EP

Premisut ORC is available in a variety of sizes and forms:

Premisut ORC			
Premisut ORC Standard	PREMISUT ORC FIBRILLAR	PREMISUT ORC KNIT	PREMISUT ORC NON-WOVEN
			
• Plant based cellulose • Fine weave of ORC • Haemostasis within 3-4 min • Completely absorbable • Bactericidal properties • Can be cut without fraying • Highly flexibility and drapeability	• Tufts of soft, light weight ORC • Activates on larger surfaces • Superior handling characteristics • Completely absorbable • Bactericidal properties • 7 layers fibrillous fabric • Useful on irregular surfaces like in Orthopaedic, Cardiovascular, Neuro and Vascular surgeries	• Denser weave of ORC for heavier bleeding • Haemostasis within 3-4 min • Bactericidal properties • Can be cut without fraying • Density suitable for suturing • Laparoscopic compatibili • Diffuse organ bleeding • Transplant procedures	• Nonwoven design increases contact surface • Needle punched with interlocking fibres • Maximized effect and superior handling • Easy use in open and MIS procedures • Ideal in Urology, Gynecology, Cardiovascular, General and Peripheral Vascular surgeries

Premisut ORC			
PREMISUT ORC TYPE	Reference	Sizes	PCS/BOX
Premisut ORC Standard	PSOO-0214	2 inch X 14 inch (5.1 cm X 35.6 cm)	12
	PSOO-0408	4 inch X 8 inch (10.2 cm X 20.3 cm)	12
	PSOO-0203	2 inch X 3 inch (5.1 cm X 7.6 cm)	12
	PSOO-502	0.5 inch X 2 inch (1.3 cm X 5.1 cm)	12
	PSOO-0102	1 inch X 2 inch (2.5 cm X 5.1 cm)	12
	PSOO-0304	3 inch X 4 inch (7.6 cm X 10.1 cm)	12
	PSOO-0205	2 inch X 5 inch (5.1 cm X 12.7 cm)	12
	PSOO-0101	1 inch X 1 inch (2.5 cm X 2.5 cm)	12
	PSOO-66	0.6 inch X 0.6 inch (1.5 cm X 1.5 cm)	12
Premisut ORC Knit	PSOK-0203	2 inch X 3 inch (5.1 cm X 7.6 cm)	12
	PSOK-0304	3 inch X 4 inch (7.6 cm X 10.2 cm)	12
	PSOK-0609	6 inch X 9 inch (15.2 cm X 22.9 cm)	12
	PSOK-0101	1 inch X 1 inch (2.5 cm X 2.5 cm)	12
	PSOK-0103	1 inch X 3 inch (2.5 cm X 7.6 cm)	12
	PSOK-0135	1 inch X 3.5 inch (2.5 cm X 8.9 cm)	12
	PSOK-0102	1 inch X 2 inch (2.5 cm X 5.1 cm)	12
	PSOK-0204	2 inch X 4 inch (5.1 cm X 10.2 cm)	12
	PSOK-0404	4 inch X 4 inch (10.2 cm X 10.2 cm)	12
Premisut ORC Fibrillar	PSOF-0102	1 inch X 2 inch (2.5 cm X 5.1 cm)	12
	PSOF-0204	2 inch X 4 inch (5.1 cm X 10.2 cm)	12
	PSOF-0404	4 inch X 4 inch (10.2 cm X 10.2 cm)	12
	PSOF-0203	2 inch X 3 inch (5.1 cm X 7.6 cm)	12
	PSOF-0304	3 inch X 4 inch (7.6 cm X 10.2 cm)	12
	PSOF-0202	2 inch X 2 inch (5.1 cm X 5.1 cm)	12
	PSOF-0101	1 inch X 1 inch (2.5 cm X 2.5 cm)	12
Premisut ORC Non-Woven	PSON-0102	1 inch X 2 inch (2.5 cm X 5.1 cm)	12
	PSON-0204	2 inch X 4 inch (5.1 cm X 10.2 cm)	12
	PSON-0404	4 inch X 4 inch (10.2 cm X 10.2 cm)	12
	PSON-0203	2 inch X 3 inch (5.1 cm X 7.6 cm)	12
	PSON-0304	3 inch X 4 inch (7.6 cm X 10.2 cm)	12
	PSON-0202	2 inch X 2 inch (5.1 cm X 5.1 cm)	12