

Premisut Spon

Premisut

Description

Premisut Spon is an Absorbable Hemostatic Gelatin Sponge, porcine origin and is a quick and effective hemostatic which stops bleeding fast. It reduces intra operative bleeding and surgery time. Premisut Spon is non-toxic, non-allergenic, non-immunogenic, and non-pyrogenic. Premisut Spon is already gamma-sterilized, and thus does not need to be re-sterilized. It is available in a variety of sizes and formats. Premisut Spon absorbs approximately 40-50 times its weight and is being completely absorbed between 3-4 weeks.

Absorption

Completely in approximately 3-4 weeks

Application areas

Abdominal, anorectal, dental surgery, ENT, urinary, gynecology, hysterectomy, orthopedic, neurosurgery, Otolaryngological, Parathyroidectomy, Vascular



Premisut Reference on CE certificate	PRODUCT NAME	PRODUCT DIMENSIONS (MM)
PMS-705010	REGULAR	70 x 50 x 10
PMS-805010	STANDARD	80 x 50 x 10
PMS-101010	DENTAL CUBES	10 x 10 x 10
PMS-1208	DENTAL TAMPON	12 x 08 Ø
PMS-202007	DENTAL	20 x 20 x 07
PMS-8030	ANAL TAMPON	80 x 30 Ø
PMS-804010	STERIFEN	80 x 40 x 10
PMS-802010	STRIP	80 x 20 x 10
PMS-802507	GYNAEC	80 x 25 x 07
PMS-801510	NASAL	80 x 15 x 07
PMS-602007	SIZE 12-7	60 x 20 x 07
PMS-100	SIZE 100	125 x 80 x 10
PMS-120	SIZE 120	120 x 100 x 10
PMS-200	SIZE 200	250 x 80 x 10
PMS-705001	SPECIAL	70 x 50 x 01
PMS-805001	SPECIAL	80 x 50 x 01
PMS-200705	FILM	200 x 70 x 0.5
PMS-8007	NASAL TAMPON	80 x 07 Ø
PMS-401010	STRIPS	40 x 10 x 10
PMS-303010	DIAL	30 x 30 x 10
PMS-01	POWDER	1 gram
PMS-02	POWDER	2 gram
PMS-03	POWDER	3 gram
PMS-05	POWDER	5 gram
PMS-01B	POWDER	1 gram
PMS-02B	POWDER	2 gram
PMS-03B	POWDER	3 gram
PMS-05B	POWDER	5 gram