

Description

Premisut NYLON suture [Polyamide(Premisut NYLON)], is a sterile synthetic monofilament non-absorbable surgical suture prepared from a copolymer of Polyamide 6 and Polyamide 6/6 (Premisut NYLON 6 and Premisut NYLON 6/6). Premisut NYLON suture is available in black and blue dyed forms.

Indications

Premisut NYLON suture is intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological tissues.

Performance

Premisut NYLON suture elicits a minimal initial inflammatory reaction in tissue which is followed by gradual encapsulation of suture by fibrous connective tissue. While Premisut NYLON suture is not absorbed, progressive hydrolysis of suture in vivo may result in gradual loss of tensile strength over a period of time.

Tensile Strength


Exceeds USP & EP requirements

Total Absorbption

None

Color

Black & Blue



Sizes

10/0 to 2 USP

Usage Areas

General surgery, skin closure, ophthalmic, oral, gastrointestinal, gynecology, plastic, microsurgery, obstetrics, neurosurgery, orthopedic.

Sterilization


Ethylene Oxide

Shelf Life

5 years

Packaging

Box of 12pcs / 36pcs



Premisut NYLON

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 16		5/0	1	Black	75	AM642HR1650	36
		4/0	1.5	Black	75	AM642HR1640	36
HR 16 PRIME NEEDLE		5/0	1	Black	75	AMP642HR1650	36
		4/0	1.5	Black	75	AMP642HR1640	36
HR 20		6/0	0.7	Black	75	AM642HR2060*	12
		5/0	1	Black	75	AM642HR2050	36
		4/0	1.5	Black	75	AM642HR2040	36
		3/0	2	Black	75	AM642HR2030	36
		2/0	3	Black	75	AM642HR2020	36
HR 20 PRIME NEEDLE		6/0	0.7	Black	75	AMP642HR2060*	12
		5/0	1	Black	75	AMP642HR2050	36
		4/0	1.5	Black	75	AMP642HR2040	36
		3/0	2	Black	75	AMP642HR2030	36
		2/0	3	Black	75	AMP642HR2020	36
HR 26		4/0	1.5	Black	75	AM642HR2640	36
		3/0	2	Black	75	AM642HR2630	36
		2/0	3	Black	75	AM642HR2620	36
		0	3.5	Black	75	AM642HR2610	36
		1	4	Black	75	AM642HR2601	36

• RIBBED
• RTT
*NO RTT

• RIBBED
• RTT
*NO RTT

• RIBBED
• RTT

Premisut NYLON

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 26 PRIME NEEDLE		4/0	1.5	Black	75	AMP642HR2640	36
		3/0	2	Black	75	AMP642HR2630	36
		2/0	3	Black	75	AMP642HR2620	36
		0	3.5	Black	75	AMP642HR2610	36
		1	4	Black	75	AMP642HR2601	36
• RIBBED • RTT							
HR 30		2/0	3	Black	75	AM642HR3020	36
		0	3.5	Black	75	AM642HR3010	36
		0	3.5	Black	150 loop	AM645HR30L10*	36
		1	4	Black	75	AM642HR3001	36
		1	4	Black	150 loop	AM645HR30L01*	36
• RIBBED • RTT *NO RTT							
HR 30 PRIME NEEDLE		2/0	3	Black	75	AMP642HR3020	36
		0	3.5	Black	75	AMP642HR3010	36
		0	3.5	Black	150 loop	AMP645HR30L10*	36
		1	4	Black	75	AMP642HR3001	36
		1	4	Black	150 loop	AMP645HR30L01*	36
• RIBBED • RTT *NO RTT							
HR 35		0	3.5	Black	150 loop	AM645HR35L10	36
		1	4	Black	150 loop	AM645HR35L01	36
• RIBBED							
HR 35 PRIME NEEDLE		0	3.5	Black	150 loop	AMP645HR35L10	36
		1	4	Black	150 loop	AMP645HR35L01	36
• RIBBED							

Premisut NYLON

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 37 <							

Premisut NYLON

● ½ Circle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 50		0	3.5	Black	150 loop	AM645HR50L10*	36
		1	4	Black	90	AM643HR5001	36
		1	4	Black	150 loop	AM645HR50L01*	36
		2	5	Black	150 loop	AM645HR50L02*	36
• RIBBED • RTT *NO RTT							
HR 50 PRIME NEEDLE		0	3.5	Black	150 loop	AMP645HR50L10*	36
		1	4	Black	90	AMP643HR5001	36
		1	4	Black	150 loop	AMP645HR50L01*	36
• RIBBED • RTT *NO RTT							

▼ 1/2 circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HS 17		4/0	1.5	Black	75	AM642HS1740	36
		3/0	2	Black	75	AM642HS1730	36
HS 17 PRIME NEEDLE		4/0	1.5	Black	75	AMP642HS1740	36
		3/0	2	Black	75	AMP642HS1730	36

Premisut NYLON

▼ 1/2 circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HS 26		3/0	2	Black	75	AM642HS2630	36
		2/0	3	Black	75	AM642HS2620	36
• RIBBED • RTT							
HS 26 PRIME NEEDLE		3/0	2	Black	75	AMP642HS2630	36
		2/0	3	Black	75	AMP642HS2620	36
• RIBBED • RTT							
HS 36		0	3.5	Black	75	AM642HS3610	36
		1	4	Black	75	AM642HS3601	36
		2	5	Black	75	AM642HS3602	36
• RIBBED • RTT							
HS 36 PRIME NEEDLE		0	3.5	Black	75	AMP642HS3610	36
		1	4	Black	75	AMP642HS3601	36
		2	5	Black	75	AMP642HS3602	36
• RIBBED • RTT							
HS 40		0	3.5	Black	75	AM642HS4010	36
		1	4	Black	75	AM642HS4001	36
		1	4	Black	100	AM647HS4001*	36
• RIBBED • RTT *NO RTT							

Premisut NYLON

▼ 1/2 circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HS 40 PRIME NEEDLE		0	3.5	Black	75	AMP642HS4010	36
		1	4	Black	75	AMP642HS4001	36
		1	4	Black	100	AMP647HS4001*	36
		• RIBBED • RTT *NO RTT					
HS 48		1	4	Black	100	AM647HS4801	36
• RIBBED							
HS 48 PRIME NEEDLE		1	4	Black	100	AMP647HS4801	36
• RIBBED							

▼ 1/2 Circle Extra Reverse Cutting – Long cutting edge (plastic/cosmetic)

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
XHS 8 PRIME NEEDLE		8/0	0.4	Black	45	AMP641XHS880	12
XHS 13 PRIME NEEDLE		5/0	1	Black	45	AMP641XHS1350	36

Premisut NYLON

▼ 1/2 Circle Extra Reverse Cutting – Long cutting edge (plastic/cosmetic) ✂

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
XHS 17 PRIME NEEDLE		5/0	1	Black	75	AMP642XHS1750	36
XHS 22 PRIME NEEDLE		3/0	2	Black	45	AMP641XHS2230	36
• RTT							
XHS 48 PRIME NEEDLE		1	4	Black	100	AMP647XHS4801	36
• RTT							

▼ 3/8 Circle Reverse Cutting ✂

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 12 PRIME NEEDLE		7/0	0.5	Black	75	AMP642DS1270	12
		6/0	0.7	Black	75	AMP642DS1260	12
		5/0	1	Black	75	AMP642DS1250	36
		4/0	1.5	Black	75	AMP642DS1240	36

Premisut NYLON

▼ 3/8 Circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 16		5/0	1	Black	45	AM641DS1650	36
		5/0	1	Black	75	AM642DS1650	36
		4/0	1.5	Black	45	AM641DS1640	36
		4/0	1.5	Black	75	AM642DS1640	36
		3/0	2	Black	75	AM642DS1630	36
DS 16 PRIME NEEDLE		6/0	0.7	Black	45	AMP641DS1660	12
		6/0	0.7	Black	75	AMP642DS1660	12
		5/0	1	Black	45	AMP641DS1650	36
		5/0	1	Black	75	AMP642DS1650	36
		4/0	1.5	Black	45	AMP641DS1640	36
		4/0	1.5	Black	75	AMP642DS1640	36
		3/0	2	Black	75	AMP642DS1630	36
DS 19 • RIBBED • RTT		5/0	1	Black	45	AM641DS1950	36
		5/0	1	Black	75	AM642DS1950	36
		4/0	1.5	Black	45	AM641DS1940	36
		4/0	1.5	Black	75	AM642DS1940	36
		3/0	2	Black	45	AM641DS1930	36
		3/0	2	Black	75	AM642DS1930	36
		2/0	3	Black	75	AM642DS1920	36
DS 19 PRIME NEEDLE • RIBBED • RTT		5/0	1	Black	45	AMP641DS1950	36
		5/0	1	Black	75	AMP642DS1950	36
		4/0	1.5	Black	45	AMP641DS1940	36
		4/0	1.5	Black	75	AMP642DS1940	36
		3/0	2	Black	45	AMP641DS1930	36
		3/0	2	Black	75	AMP642DS1930	36
		2/0	3	Black	75	AMP642DS1920	36
DS 22 • RIBBED • RTT		4/0	1.5	Black	75	AM642DS2240	36

Premisut NYLON

▼ 3/8 Circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 22 PRIME NEEDLE		4/0	1.5	Black	75	AMP642DS2240	36
• RIBBED • RTT							
DS 24		5/0	1	Black	45	AM641DS2450	36
		5/0	1	Black	75	AM642DS2450	36
		4/0	1.5	Black	45	AM641DS2440	36
		4/0	1.5	Black	75	AM642DS2440	36
		3/0	2	Black	45	AM641DS2430	36
		3/0	2	Black	75	AM642DS2430	36
		2/0	3	Black	45	AM641DS2420	36
		2/0	3	Black	75	AM642DS2420	36
		0	3.5	Black	75	AM642DS2410	36
		1	4	Black	75	AM642DS2401	36
		2	5	Black	75	AM642DS2402	36
		2	5	Black	90	AM643DS2402*	36
• RIBBED • RTT *NO RTT							
DS 24 PRIME NEEDLE		5/0	1	Black	45	AMP641DS2450	36
		5/0	1	Black	75	AMP642DS2450	36
		4/0	1.5	Black	45	AMP641DS2440	36
		4/0	1.5	Black	75	AMP642DS2440	36
		3/0	2	Black	45	AMP641DS2430	36
		3/0	2	Black	75	AMP642DS2430	36
		2/0	3	Black	45	AMP641DS2420	36
		2/0	3	Black	75	AMP642DS2420	36
		0	3.5	Black	75	AMP642DS2410	36
		1	4	Black	75	AMP642DS2401	36
		2	5	Black	75	AMP642DS2402	36
		2	5	Black	90	AMP643DS2402*	36
• RIBBED • RTT *NO RTT							

Premisut NYLON

▼ 3/8 Circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 26		4/0	1.5	Black	45	AM641DS2640	36
		3/0	2	Black	45	AM641DS2630	36
		3/0	2	Black	75	AM642DS2630	36
		2/0	3	Black	45	AM641DS2620	36
		2/0	3	Black	75	AM642DS2620	36
• RIBBED • RTT							
DS 26 PRIME NEEDLE		4/0	1.5	Black	45	AMP641DS2640	36
		3/0	2	Black	45	AMP641DS2630	36
		3/0	2	Black	75	AMP642DS2630	36
		2/0	3	Black	45	AMP641DS2620	36
		2/0	3	Black	75	AMP642DS2620	36
• RIBBED • RTT							
DS 30		4/0	1.5	Black	75	AM642DS3040	36
		3/0	2	Black	45	AM641DS3030	36
		3/0	2	Black	75	AM642DS3030	36
		2/0	3	Black	75	AM642DS3020	36
		0	3.5	Black	75	AM642DS3010	36
		1	4	Black	75	AM642DS3001	36
• RIBBED • RTT							
DS 30 PRIME NEEDLE		4/0	1.5	Black	75	AMP642DS3040	36
		3/0	2	Black	45	AMP641DS3030	36
		3/0	2	Black	75	AMP642DS3030	36
		2/0	3	Black	75	AMP642DS3020	36
		0	3.5	Black	75	AMP642DS3010	36
		1	4	Black	75	AMP642DS3001	36
• RIBBED • RTT							
DS 35		3/0	2	Black	75	AM642DS3530	36
		2/0	3	Black	75	AM642DS3520	36
		0	3.5	Black	75	AM642DS3510	36
• RIBBED • RTT							

Premisut NYLON

▼ 3/8 Circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 35 PRIME NEEDLE		3/0	2	Black	75	AMP642DS3530	36
		2/0	3	Black	75	AMP642DS3520	36
		0	3.5	Black	75	AMP642DS3510	36
DS 40 PRIME NEEDLE		2/0	3	Black	90	AM643DS4020	36
		0	3.5	Black	90	AM643DS4010	36
		1	4	Black	90	AM643DS4001	36
		2	5	Black	90	AM643DS4002*	36
DS 40 PRIME NEEDLE		2/0	3	Black	90	AMP643DS4020	36
		0	3.5	Black	90	AMP643DS4010	36
		1	4	Black	90	AMP643DS4001	36
		2	5	Black	90	AMP643DS4002*	36
DS 45 PRIME NEEDLE		2/0	3	Black	90	AM643DS4520	36
		0	3.5	Black	90	AM643DS4510	36
		1	4	Black	90	AM643DS4501	36
DS 45 PRIME NEEDLE		2/0	3	Black	90	AMP643DS4520	36
		0	3.5	Black	90	AMP643DS4510	36
		1	4	Black	90	AMP643DS4501	36

Premisut NYLON

▼ 3/8 Circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 50		0	3.5	Black	90	AM643DS5010	36
		1	4	Black	90	AM643DS5001	36
• RIBBED							
DS 50 PRIME NEEDLE		0	3.5	Black	90	AMP643DS5010	36
		1	4	Black	90	AMP643DS5001	36
• RIBBED							
DS 60		2/0	3	Black	75	AM642DS6020	36
		2/0	3	Black	100	AM647DS6020	36
		0	3.5	Black	75	AM642DS6010	36
• RIBBED							
DS 60 PRIME NEEDLE		2/0	3	Black	75	AMP642DS6020	36
		2/0	3	Black	100	AMP647DS6020	36
		0	3.5	Black	75	AMP642DS6010	36
• RIBBED							
DS 90		3/0	2	Black	90	AM643DS9030	36
		2/0	3	Black	90	AM643DS9020	36
		0	3.5	Black	90	AM643DS9010	36
		1	4	Black	90	AM643DS9001	36
		2	5	Black	90	AM643DS9002	36
• RIBBED							

Premisut NYLON

3/8 Circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 90 PRIME NEEDLE		3/0	2	Black	90	AMP643DS9030	36
		2/0	3	Black	90	AMP643DS9020	36
		0	3.5	Black	90	AMP643DS9010	36
		1	4	Black	90	AMP643DS9001	36
		2	5	Black	90	AMP643DS9002	36
• RIBBED							

3/8 Circle Reverse Cutting – Long cutting edge (Extra Cutting/cosmetic)

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
XDS 8 PRIME NEEDLE		7/0	0.5	Black	45	AMP641XDS870	12
XDS 11 PRIME NEEDLE		7/0	0.5	Black	45	AMP641XDS1170	12
		7/0	0.5	Black	75	AMP642XDS1170	12
		6/0	0.7	Black	45	AMP641XDS1160	12
		6/0	0.7	Black	75	AMP642XDS1160	12
XDS 13 PRIME NEEDLE		6/0	0.7	Black	45	AMP641XDS1360	12
		6/0	0.7	Black	75	AMP642XDS1360	12
		5/0	1	Black	45	AMP641XDS1350	36
		5/0	1	Black	75	AMP642XDS1350	36
		4/0	1.5	Black	45	AMP641XDS1340	36
		4/0	1.5	Black	75	AMP642XDS1340	36

Premisut NYLON

▼ **3/8 Circle Reverse Cutting – Long cutting edge (Extra Cutting/cosmetic)** 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
XDS 16 PRIME NEEDLE		6/0	0.7	Black	45	AMP641XDS1660	12
		6/0	0.7	Black	75	AMP642XDS1660	12
		5/0	1	Black	45	AMP641XDS1650	36
		5/0	1	Black	75	AMP642XDS1660	36
		4/0	1.5	Black	45	AMP641XDS1640	36
		4/0	1.5	Black	75	AMP642XDS1640	36
XDS 19 PRIME NEEDLE		5/0	1	Black	45	AMP641XDS1950	36
		5/0	1	Black	75	AMP642XDS1950	36
		4/0	1.5	Black	45	AMP641XDS1940	36
		3/0	2	Black	45	AMP641XDS1930	36
		3/0	2	Black	75	AMP642XDS1930	36
• RTT							
XDS 24 PRIME NEEDLE		4/0	1.5	Black	45	AMP641XDS2440	36
		4/0	1.5	Black	75	AMP642XDS2440	36
		3/0	2	Black	45	AMP641XDS2430	36
		3/0	2	Black	75	AMP642XDS2430	36
• RTT							
XDS 26 PRIME NEEDLE		3/0	2	Black	45	AMP641XDS2630	36
		3/0	2	Black	75	AMP642XDS2630	36
• RTT							
XDS 30 PRIME NEEDLE		3/0	2	Black	75	AMP642XDS3030	36
• RTT							

Premisut NYLON

▲ 3/8 Circle Conventional Long cutting edge (Extra Cutting/cosmetic) ▲

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
XNS 13 PRIME NEEDLE		6/0	0.7	Black	45	AMP641XNS1360	12
XNS 16 PRIME NEEDLE		6/0	0.7	Black	45	AMP641XNS1660	12
		6/0	0.7	Black	75	AMP642XNS1660	12
		5/0	1	Black	45	AMP641XNS1650	36
		5/0	1	Black	75	AMP642XNS1650	36
XNS 19 PRIME NEEDLE		4/0	1.5	Black	45	AMP641XNS1940	36
		4/0	1.5	Black	75	AMP642XNS1940	36
• RTT							
XNS 26 PRIME NEEDLE		6/0	0.7	Black	45	AMP641XNS2660*	12
		5/0	1	Black	45	AMP641XNS2650	36
• RTT *NO RTT							

● 3/8 Circle Round bodied TAPERPOINT ▲

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DR 3.8 PRIME NEEDLE		10/0	0.2	Black	15	AMP649DR3.8100	12
		9/0	0.3	Black	15	AMP649DR3.890	12

Premisut NYLON

● 3/8 Circle Round bodied TAPERPOINT 


NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DR 4.7 PRIME NEEDLE		9/0	0.3	Black	15	AMP649DR4.790	12
		8/0	0.4	Black	15	AMP649DR4.780	12
DR 5.1 PRIME NEEDLE		10/0	0.2	Black	15	AMP649DR5.1100	12
		9/0	0.3	Black	15	AMP649DR5.190	12
DR 6.5 PRIME NEEDLE		10/0	0.2	Black	15	AMP649DR6.5100	12
		9/0	0.3	Black	15	AMP649DR6.590	12
		8/0	0.4	Black	15	AMP649DR6.580	12
2 x DR 6.5 PRIME NEEDLE		10/0	0.2	Black	30	AMP644DR6.5100-2	12
DR 8 PRIME NEEDLE		8/0	0.4	Black	45	AMP641DR880	12

Premisut NYLON

3/8 Circle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DR 9.3 PRIME NEEDLE		7/0	0.5	Black	45	AMP641DR9.370	12
DR 13 PRIME NEEDLE		6/0	0.7	Black	75	AMP642DR1360	12

Straight needle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2 x GR 65 PRIME NEEDLE		2/0	3	Black	75	AMP642GR6520-2	36

Straight Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
GS 50 • RIBBED		2/0	3	Black	75	AM642GS5020	36

Premisut NYLON

▲ Straight Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
GS 50 PRIME NEEDLE		2/0	3	Black	75	AMP642GS5020	36
• RIBBED							
GS 60		3/0	2	Black	75	AM642GS6030	36
		2/0	3	Black	75	AM642GS6020	36
• RIBBED							
GS 60 PRIME NEEDLE		3/0	2	Black	75	AMP642GS6030	36
		2/0	3	Black	75	AMP642GS6020	36
• RIBBED							

▼ 3/8 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DSP 6 PRIME NEEDLE		10/0	0.2	Black	15	AMP649DSP6100	12

Premisut NYLON

3/8 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2 x DSP 6.5 PRIME NEEDLE		10/0	0.2	Black	30	AMP644DSP6.5100-2	12
		9/0	0.3	Black	30	AMP644DSP6.590-2	12

1/2 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2 x HSP 5.5 PRIME NEEDLE		10/0	0.2	Black	30	AMP644HSP5.5100-2	12
2 x HSP 7 PRIME NEEDLE		8/0	0.4	Black	30	AMP644HSP780-2	12

1/4 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2 x VSP 6.5 PRIME NEEDLE		8/0	0.4	Black	30	AMP644VSP6.580-2	12

Premisut NYLON

1/4 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
VSP 8 PRIME NEEDLE		5/0	1	Black	30	AMP644VSP850	36

Without needle

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
	without needle 3 x 45cm	4/0	1.5	Black	3 x 45	AM64-3.4540	36
		3/0	2	Black	3 x 45	AM64-3.4530	36
		2/0	3	Black	3 x 45	AM64-3.4520	36
		0	3.5	Black	3 x 45	AM64-3.4510	36
	without needle 10 x 45cm	3/0	2	Black	10 x 45	AM64-10.4530	36
12 x 45cm		4/0	1.5	Black	12 x 45	AM64-12.4540	36
		3/0	2	Black	12 x 45	AM64-12.4530	36
		0	3.5	Black	12 x 45	AM64-12.4510	36
	without needle 2 x 70cm	2/0	3	Black	2 x 70	AM64-2.7020	36
	without needle 10 x 70cm	2/0	3	Black	10 x 75	AM64-10.7520	36