

Premisut PGLA-910

Premisut

Description

Premisut PGLA-910 suture (Polyglactin 910), is a sterile synthetic multifilament braided absorbable surgical suture prepared from a co-polymer of 90% Glycolide and 10% L-lactide. Premisut PGLA-910 suture is dyed with D&C Violet No 2 (conforming to U.S.Code of Federal regulations). Premisut PGLA- 910 suture is available in undyed form also. For added lubrication, smooth passage through tissue and easy knot tie down, Premisut PGLA-910 suture is coated with a mixture composed of Polyglactin 370 co-polymer (30% Glycolide and 70% Lactide) and Calcium Stearate.

Indications

Premisut PGLA-910 suture is intended for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological issues

Performance

Premisut PGLA-910 suture after implantation elicits a minimal initial inflammatory reaction in tissues and is eventually replaced with an ingrowth of fibrous connective tissue. Premisut PGLA-910 suture being absorbable in nature gets absorbed in the tissues by means of hydrolysis, during which the Polyglactin 910 co-polymer is degraded into Glycolic acid and Lactic acid which is eventually absorbed and metabolized in the body. This process of absorption initially leads to progressive loss of tensile strength, followed by a loss of mass.

Days post implantation

14 days
21 days
28 days

Tensile strength

75%
49%
27%

Total Absorption

~60~70 days

Color

Violet & Undyed



Sizes

8/0 to 2 USP

Usage Areas

General surgery, general soft tissue approximation and/or ligation, gastrointestinal, urology, orthopedic, ophthalmic, gynecology, skin closure, plastic surgery, obstetrics, peripheral nerve anastomosis, microsurgery.

Sterilization

Ethylene Oxide

Shelf Life

5 years

Packaging

Box of 12pcs / 36pcs



Premisut PGLA - 910

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 10 PRIME NEEDLE		5/0	1	Violet	70	AMP1L18HR1050	36
HR 13		6/0	0.7	Violet	75	AM1L12HR1360	12
		5/0	1	Violet	75	AM1L12HR1350	36
		4/0	1.5	Violet	70	AM1L18HR1340	36
		3/0	2	Violet	70	AM1L18HR1330	36
		2/0	3	Violet	70	AM1L18HR1320	36
		0	3.5	Violet	70	AM1L18HR1310	36
HR 13 PRIME NEEDLE		6/0	0.7	Violet	75	AMP1L12HR1360	12
		5/0	1	Violet	75	AMP1L12HR1350	36
		4/0	1.5	Violet	70	AMP1L18HR1340	36
		3/0	2	Violet	70	AMP1L18HR1330	36
		2/0	3	Violet	70	AMP1L18HR1320	36
		0	3.5	Violet	70	AMP1L18HR1310	36
HR 17		6/0	0.7	Violet	70	AM1L18HR1760	12
		5/0	1	Violet	75	AM1L12HR1750	36
		5/0	1	Undyed	70	AM1L28HR1750	36
		4/0	1.5	Violet	75	AM1L12HR1740	36
		4/0	1.5	Undyed	70	AM1L28HR1740	36
		3/0	2	Violet	75	AM1L12HR1730	36
		3/0	2	Undyed	70	AM1L28HR1730	36
		2/0	3	Violet	75	AM1L12HR1720	36
HR 17 PRIME NEEDLE		6/0	0.7	Violet	70	AMP1L18HR1760	12
		5/0	1	Violet	75	AMP1L12HR1750	36
		5/0	1	Undyed	70	AMP1L28HR1750	36
		4/0	1.5	Violet	75	AMP1L12HR1740	36
		4/0	1.5	Undyed	70	AMP1L28HR1740	36
		3/0	2	Violet	75	AMP1L12HR1730	36
		3/0	2	Undyed	70	AMP1L28HR1730	36
		2/0	3	Violet	75	AMP1L12HR1720	36

Premisut PGLA - 910

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2xHR 17		3/0	2	Violet	70	AM1L18HR1730-2	36
2xHR 17 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18HR1730-2	36
HR 17 VisiBlack PRIME NEEDLE		5/0	1	Violet	70	AMP1L18HR1750VB	36
		4/0	1.5	Violet	70	AMP1L18HR1740VB	36
HR 20 • RIBBED • RTT		4/0	1.5	Violet	75	AM1L12HR2040	36
		3/0	2	Violet	75	AM1L12HR2030	36
		2/0	3	Violet	75	AM1L12HR2020	36
HR 20 PRIME NEEDLE • RIBBED • RTT		4/0	1.5	Violet	75	AMP1L12HR2040	36
		3/0	2	Violet	75	AMP1L12HR2030	36
		2/0	3	Violet	75	AMP1L12HR2020	36

Premisut PGLA - 910

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM			USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 22 • RIBBED • RTT			4/0	1.5	Violet	75	AM1L12HR2240	36
			4/0	1.5	Undyed	70	AM1L28HR2240	36
			3/0	2	Violet	75	AM1L12HR2230	36
			3/0	2	Undyed	70	AM1L28HR2230	36
			2/0	3	Violet	75	AM1L12HR2220	36
			2/0	3	Undyed	70	AM1L28HR2220	36
HR 22 PRIME NEEDLE • RIBBED • RTT			4/0	1.5	Violet	75	AMP1L12HR2240	36
			4/0	1.5	Undyed	70	AMP1L28HR2240	36
			3/0	2	Violet	75	AMP1L12HR2230	36
			3/0	2	Undyed	70	AMP1L28HR2230	36
			2/0	3	Violet	75	AMP1L12HR2220	36
			2/0	3	Undyed	70	AMP1L28HR2220	36
HR 22 VisiBlack PRIME NEEDLE • RIBBED • RTT			4/0	1.5	Violet	70	AMP1L18HR2240VB	36
			3/0	2	Violet	70	AMP1L18HR2230VB	36
			2/0	3	Violet	70	AMP1L18HR2220VB	36
HR 26 								

Premisut PGLA - 910

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 26 PRIME NEEDLE		4/0	1.5	Violet	70	AMP1L18HR2640	36
		3/0	2	Violet	75	AMP1L12HR2630	36
		3/0	2	Undyed	70	AMP1L28HR2630	36
		3/0	2	Violet	90	AMP1L13HR2630	36
		2/0	3	Violet	75	AMP1L12HR2620	36
		2/0	3	Undyed	70	AMP1L28HR2620	36
		2/0	3	Violet	90	AMP1L13HR2620	36
		0	3.5	Violet	75	AMP1L12HR2610	36
		0	3.5	Undyed	70	AMP1L28HR2610	36
		1	4	Violet	75	AMP1L12HR2601	36
		2	5	Violet	75	AMP1L12HR2602	36
• RIBBED • RTT							
HR 26 VisiBlack PRIME NEEDLE		4/0	1.5	Violet	70	AMP1L18HR2640VB	36
		3/0	2	Violet	70	AMP1L18HR2630VB	36
		3/0	2	Undyed	70	AMP1L28HR2630VB	36
		2/0	3	Violet	70	AMP1L18HR2620VB	36
• RIBBED • RTT							
HR 30		3/0	2	Violet	75	AM1L12HR3030	36
		2/0	3	Violet	75	AM1L12HR3020	36
		0	3.5	Violet	75	AM1L12HR3010	36
		1	4	Violet	75	AM1L12HR3001	36
• RIBBED • RTT							
HR 30 PRIME NEEDLE		3/0	2	Violet	75	AMP1L12HR3030	36
		2/0	3	Violet	75	AMP1L12HR3020	36
		0	3.5	Violet	75	AMP1L12HR3010	36
		1	4	Violet	75	AMP1L12HR3001	36
• RIBBED • RTT							

Premisut PGLA - 910

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 30 heavy		0	3.5	Violet	75	AM1L12HR3010H	36
		1	4	Violet	75	AM1L12HR3001H	36
		• RIBBED • RTT • HEAVY NEEDLE					
HR 30 heavy PRIME NEEDLE		0	3.5	Violet	75	AMP1L12HR3010H	36
		1	4	Violet	75	AMP1L12HR3001H	36
		• RIBBED • RTT • HEAVY NEEDLE					
HR 36		4/0	1.5	Violet	90	AM1L13HR3640	36
		3/0	2	Violet	75	AM1L12HR3630	36
		3/0	2	Undyed	70	AM1L28HR3630	36
		3/0	2	Violet	90	AM1L13HR3630	36
		2/0	3	Violet	75	AM1L12HR3620	36
		2/0	3	Undyed	70	AM1L28HR3620	36
		2/0	3	Violet	90	AM1L13HR3620	36
		2/0	3	Undyed	90	AM1L23HR3620	36
		0	3.5	Violet	70	AM1L18HR3610	36
		0	3.5	Violet	90	AM1L13HR3610	36
		0	3.5	Undyed	90	AM1L23HR3610	36
		1	4	Violet	70	AM1L18HR3601	36
		1	4	Violet	90	AM1L13HR3601	36
		1	4	Undyed	70	AM1L28HR3601	36
		2	5	Violet	90	AM1L13HR3602*	36
		2	5	Undyed	90	AM1L23HR3602*	36
• RIBBED • RTT *NO RTT							

Premisut PGLA - 910

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 36 PRIME NEEDLE		4/0	1.5	Violet	90	AMP1L13HR3640	36
		3/0	2	Violet	75	AMP1L12HR3630	36
		3/0	2	Undyed	70	AMP1L28HR3630	36
		3/0	2	Violet	90	AMP1L13HR3630	36
		2/0	3	Violet	75	AMP1L12HR3620	36
		2/0	3	Undyed	70	AMP1L28HR3620	36
		2/0	3	Violet	90	AMP1L13HR3620	36
		2/0	3	Undyed	90	AMP1L23HR3620	36
		0	3.5	Violet	70	AMP1L18HR3610	36
		0	3.5	Violet	90	AMP1L13HR3610	36
		0	3.5	Undyed	90	AMP1L23HR3610	36
		1	4	Violet	70	AMP1L18HR3601	36
		1	4	Violet	90	AMP1L13HR3601	36
		1	4	Undyed	70	AMP1L28HR3601	36
		2	5	Violet	90	AMP1L13HR3602*	36
		2	5	Undyed	90	AMP1L23HR3602*	36
• RIBBED • RTT * NO RTT							
HR 40		3/0	2	Violet	70	AM1L18HR4030	36
		2/0	3	Violet	75	AM1L12HR4020	36
		0	3.5	Violet	70	AM1L18HR4010	36
		0	3.5	Violet	90	AM1L13HR4010	36
		0	3.5	Violet	150 loop	AM1L15HR40L10*	36
		1	4	Violet	70	AM1L18HR4001	36
		1	4	Violet	90	AM1L13HR4001	36
		1	4	Violet	100	AM1L17HR4001	36
		2	5	Violet	70	AM1L18HR4002	36
		2	5	Violet	90	AM1L13HR4002*	36
		2	5	Undyed	90	AM1L23HR4002*	36
• RIBBED • RTT * NO RTT							
HR 40 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18HR4030	36
		2/0	3	Violet	75	AMP1L12HR4020	36
		0	3.5	Violet	70	AMP1L18HR4010	36
		0	3.5	Violet	90	AMP1L13HR4010	36
		0	3.5	Violet	150 loop	AMP1L15HR40L10*	36
		1	4	Violet	70	AMP1L18HR4001	36
		1	4	Violet	90	AMP1L13HR4001	36
		2	5	Violet	70	AMP1L18HR4002	36
		2	5	Violet	90	AMP1L13HR4002*	36
		2	5	Undyed	90	AMP1L23HR4002*	36
• RIBBED • RTT * NO RTT							

Premisut PGLA - 910

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 40 heavy		1	4	Violet	75	AM1L12HR4001H	36
		2	5	Violet	75	AM1L12HR4002H	36
• RIBBED • RTT • HEAVY NEEDLE							
HR 40 heavy PRIME NEEDLE		1	4	Violet	75	AMP1L12HR4001H	36
		2	5	Violet	75	AMP1L12HR4002H	36
• RIBBED • RTT • HEAVY NEEDLE							
HR 45		2/0	3	Violet	75	AM1L12HR4520	36
		0	3.5	Violet	75	AM1L12HR4510	36
		1	4	Violet	70	AM1L18HR4501	36
		1	4	Violet	90	AM1L13HR4501	36
		2	5	Violet	70	AM1L18HR4502	36
		2	5	Violet	90	AM1L13HR4502*	36
• RIBBED • RTT *NO RTT							
HR 45 PRIME NEEDLE		2/0	3	Violet	75	AMP1L12HR4520	36
		0	3.5	Violet	75	AMP1L12HR4510	36
		1	4	Violet	70	AMP1L18HR4501	36
		1	4	Violet	90	AMP1L13HR4501	36
		2	5	Violet	70	AMP1L18HR4502	36
		2	5	Violet	90	AMP1L13HR4502*	36
• RIBBED • RTT *NO RTT							
HR 45 Heavy		1	4	Violet	75	AM1L12HR4501H	36
• RIBBED • RTT • HEAVY NEEDLE							

Premisut PGLA - 910

● ½ Circle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 45 Heavy PRIME NEEDLE		1	4	Violet	75	AMP1L12HR4501H	36
• RIBBED • RTT • HEAVY NEEDLE							
HR 48		3/0	2	Violet	70	AM1L18HR4830	36
		2/0	3	Violet	75	AM1L12HR4820	36
		0	3.5	Violet	75	AM1L12HR4810	36
		0	3.5	Violet	90	AM1L13HR4810	36
		0	3.5	Violet	150 loop	AM1L15HR48L10	36
		1	4	Violet	70	AM1L18HR4801	36
		1	4	Violet	90	AM1L13HR4801	36
		2	5	Violet	70	AM1L18HR4802	36
		2	5	Violet	90	AM1L13HR4802	36
		2	5	Undyed	90	AM1L23HR4802	36
• RIBBED							
HR 48 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18HR4830	36
		2/0	3	Violet	75	AMP1L12HR4820	36
		0	3.5	Violet	75	AMP1L12HR4810	36
		0	3.5	Violet	90	AMP1L13HR4810	36
		0	3.5	Violet	150 loop	AMP1L15HR48L10	36
		1	4	Violet	70	AMP1L18HR4801	36
		1	4	Violet	90	AMP1L13HR4801	36
		2	5	Violet	70	AMP1L18HR4802	36
		2	5	Violet	90	AMP1L13HR4802	36
		2	5	Undyed	90	AMP1L23HR4802	36
• RIBBED							

● ½ Circle Round bodied TAPERPOINT Ultra Plus

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRP 6.5 PRIME NEEDLE		8/0	0.4	Violet	45	AMP1L11HRP6.580	12
		7/0	0.5	Violet	45	AMP1L11HRP6.570	12

Premisut PGLA - 910

 ½ Circle Round bodied TAPERPOINT Ultra Plus 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRP 10 PRIME NEEDLE		5/0	1	Violet	70	AMP1L18HRP1050	12
HRP 13 PRIME NEEDLE		6/0	0.7	Violet	70	AMP1L18HRP1360	12
		6/0	0.7	Undyed	45	AMP1L21HRP1360	12
		5/0	1	Violet	70	AMP1L18HRP1350	36
		4/0	1.5	Violet	70	AMP1L18HRP1340	36
		3/0	2	Violet	70	AMP1L18HRP1330	36
		0	3.5	Violet	70	AMP1L18HRP1310	36
HRP 17 PRIME NEEDLE		6/0	0.7	Violet	70	AMP1L18HRP1760	12
		5/0	1	Violet	70	AMP1L18HRP1750	36
		5/0	1	Undyed	70	AMP1L28HRP1750	36
		4/0	1.5	Violet	70	AMP1L18HRP1740	36
		4/0	1.5	Undyed	70	AMP1L28HRP1740	36
		3/0	2	Violet	70	AMP1L18HRP1730	36
		3/0	2	Undyed	70	AMP1L28HRP1730	36
		2/0	3	Violet	70	AMP1L18HRP1720	36
		0	3.5	Violet	70	AMP1L18HRP1710	36
		1	4	Violet	70	AMP1L18HRP1701	36
2 X HRP 17 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18HRP1730-2	36

Premisut PGLA - 910

 ½ Circle Round bodied TAPERPOINT Ultra Plus 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRP 22 PRIME NEEDLE		4/0	1.5	Violet	70	AMP1L18HRP2240	36
		4/0	1.5	Undyed	70	AMP1L28HRP2240	36
		3/0	2	Violet	70	AMP1L18HRP2230	36
		3/0	2	Undyed	70	AMP1L28HRP2230	36
		2/0	3	Violet	70	AMP1L18HRP2220	36
		2/0	3	Undyed	70	AMP1L28HRP2220	36
		0	3.5	Violet	70	AMP1L18HRP2210	36
• RIBBED • RTT							
HRP 26 PRIME NEEDLE		4/0	1.5	Violet	70	AMP1L18HRP2640	36
		3/0	2	Violet	70	AMP1L18HRP2630	36
		3/0	2	Undyed	70	AMP1L28HRP2630	36
		3/0	2	Violet	90	AMP1L13HRP2630	36
		2/0	3	Violet	70	AMP1L18HRP2620	36
		2/0	3	Undyed	70	AMP1L28HRP2620	36
		2/0	3	Violet	90	AMP1L13HRP2620	36
		0	3.5	Violet	70	AMP1L18HRP2610	36
		0	3.5	Undyed	70	AMP1L28HRP2610	36
		1	4	Violet	70	AMP1L18HRP2601	36
		2	5	Violet	70	AMP1L18HRP2602	36
• RIBBED • RTT							
HRP 30 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18HRP3030	36
		2/0	3	Violet	70	AMP1L18HRP3020	36
		0	3.5	Violet	70	AMP1L18HRP3010	36
		1	4	Violet	75	AMP1L12HRP3001	36
• RIBBED • RTT							
HRP 36 PRIME NEEDLE		4/0	1.5	Violet	90	AMP1L13HRP3640	36
		3/0	2	Violet	70	AMP1L18HRP3630	36
		3/0	2	Violet	90	AMP1L13HRP3630	36
		2/0	3	Violet	70	AMP1L18HRP3620	36
		2/0	3	Undyed	70	AMP1L28HRP3620	36
		2/0	3	Violet	90	AMP1L13HRP3620	36
		0	3.5	Violet	70	AMP1L18HRP3610	36
		0	3.5	Violet	90	AMP1L13HRP3610	36
		1	4	Violet	70	AMP1L18HRP3601	36
		1	4	Undyed	70	AMP1L28HRP3601	36
		2	5	Violet	90	AMP1L13HRP3602*	36
• RIBBED • RTT *NO RTT							

Premisut PGLA - 910

 ½ Circle Round bodied TAPERPOINT Ultra Plus 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRP 40 PRIME NEEDLE • RIBBED • RTT *NO RTT		3/0	2	Violet	70	AMP1L18HRP4030	36
		2/0	3	Violet	70	AMP1L18HRP4020	36
		0	3.5	Violet	70	AMP1L18HRP4010	36
		0	3.5	Violet	90	AMP1L13HRP4010	36
		1	4	Violet	70	AMP1L18HRP4001	36
		1	4	Violet	90	AMP1L13HRP4001	36
		2	5	Violet	90	AMP1L13HRP4002*	36
		2	5	Undyed	90	AMP1L23HRP4002*	36
HRP 45 PRIME NEEDLE • RIBBED • RTT		1	4	Violet	75	AMP1L12HRP4501	36
HRP 48 PRIME NEEDLE • RIBBED		3/0	2	Violet	70	AMP1L18HRP4830	36
		2/0	3	Violet	70	AMP1L18HRP4820	36
		0	3.5	Violet	70	AMP1L18HRP4810	36
		0	3.5	Violet	90	AMP1L13HRP4810	36
		0	3.5	Undyed	90	AMP1L23HRP4810	36
		1	4	Violet	70	AMP1L18HRP4801	36
		1	4	Violet	90	AMP1L13HRP4801	36
		1	4	Undyed	90	AMP1L23HRP4801	36
		2	5	Violet	90	AMP1L13HRP4802	36
		2	5	Undyed	90	AMP1L23HRP4802	36

Premisut PGLA - 910

 **1/2 Circle Round bodied Blunt Point** 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRN 30		3/0	2	Violet	75	AM1L12HRN3030	36
		2/0	3	Violet	75	AM1L12HRN3020	36
• RIBBED • RTT							
HRN 30 PRIME NEEDLE		3/0	2	Violet	75	AMP1L12HRN3030	36
		2/0	3	Violet	75	AMP1L12HRN3020	36
• RIBBED • RTT							
HRN 36 PRIME NEEDLE		2/0	3	Violet	75	AMP1L12HRN3620	36
		0	3.5	Violet	75	AMP1L12HRN3610	36
• RIBBED • RTT							
HRN 36 heavy		2/0	3	Violet	75	AM1L12HRN3620H	36
		0	3.5	Violet	90	AM1L13HRN3610H	36
		1	4	Violet	90	AM1L13HRN3601H	36
• RIBBED • RTT • HEAVY NEEDLE							
HRN 36 heavy PRIME NEEDLE		2/0	3	Violet	75	AMP1L12HRN3620H	36
		0	3.5	Violet	90	AMP1L13HRN3610H	36
		1	4	Violet	90	AMP1L13HRN3601H	36
• RIBBED • RTT • HEAVY NEEDLE							

Premisut PGLA - 910

○ 1/2 Circle Round bodied Blunt Point 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRN 40		0	3.5	Violet	75	AM1L12HRN4010	36
• RIBBED • RTT							
HRN 40 PRIME NEEDLE		0	3.5	Violet	75	AMP1L12HRN4010	36
• RIBBED • RTT							
HRN 40 heavy		0	3.5	Violet	90	AM1L13HRN4010H	36
		1	4	Violet	90	AM1L13HRN4001H	36
• RIBBED • RTT • HEAVY NEEDLE							
HRN 40 heavy PRIME NEEDLE		0	3.5	Violet	90	AMP1L13HRN4010H	36
		1	4	Violet	90	AMP1L13HRN4001H	36
• RIBBED • RTT • HEAVY NEEDLE							
HRN 48		0	3.5	Violet	90	AM1L13HRN4810	36
		1	4	Violet	90	AM1L13HRN4801	36
• RIBBED							

Premisut PGLA - 910

1/2 Circle Round bodied Blunt Point

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRN 48 PRIME NEEDLE		0	3.5	Violet	90	AMP1L13HRN4810	36
		1	4	Violet	90	AMP1L13HRN4801	36
• RIBBED							
HRN 65		1	4	Violet	100	AM1L17HRN6501	36
• RIBBED							
HRN 65 PRIME NEEDLE		1	4	Violet	100	AMP1L17HRN6501	36
• RIBBED							

1/2 circle round bodied TAPERCUT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRT 17		4/0	1.5	Violet	70	AM1L18HRT1740	36
		4/0	1.5	Undyed	70	AM1L28HRT1740	36
		3/0	2	Violet	70	AM1L18HRT1730	36
HRT 17 PRIME NEEDLE		4/0	1.5	Violet	70	AMP1L18HRT1740	36
		4/0	1.5	Undyed	70	AMP1L28HRT1740	36
		3/0	2	Violet	70	AMP1L18HRT1730	36

Premisut PGLA - 910

 **1/2 circle round bodied TAPERCUT** 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRT 26		3/0	2	Violet	70	AM1L18HRT2630	36
		2/0	3	Violet	70	AM1L18HRT2620	36
• RIBBED • RTT							
HRT 26 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18HRT2630	36
		2/0	3	Violet	70	AMP1L18HRT2620	36
• RIBBED • RTT							
HRT 36		3/0	2	Violet	90	AM1L13HRT3630	36
		2/0	3	Violet	90	AM1L13HRT3620	36
		0	3.5	Violet	90	AM1L13HRT3610	36
		1	4	Violet	90	AM1L13HRT3601	36
		1	4	Undyed	90	AM1L23HRT3601	36
• RIBBED • RTT							
HRT 36 PRIME NEEDLE		3/0	2	Violet	90	AMP1L13HRT3630	36
		2/0	3	Violet	90	AMP1L13HRT3620	36
		0	3.5	Violet	90	AMP1L13HRT3610	36
		1	4	Violet	90	AMP1L13HRT3601	36
		1	4	Undyed	90	AMP1L23HRT3601	36
• RIBBED • RTT							
HRT 40		2/0	3	Violet	90	AM1L13HRT4020	36
		0	3.5	Violet	90	AM1L13HRT4010	36
• RIBBED • RTT							

Premisut PGLA - 910

 1/2 circle round bodied TAPERCUT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRT 40 PRIME NEEDLE		2/0	3	Violet	90	AMP1L13HRT4020	36
		0	3.5	Violet	90	AMP1L13HRT4010	36
• RIBBED • RTT							
2xHRT 40		1	4	Violet	70	AM1L18HRT4001-2	36
• RIBBED							
2xHRT 40 PRIME NEEDLE		1	4	Violet	70	AMP1L18HRT4001-2	36
• RIBBED							
HRT 48		2/0	3	Violet	90	AM1L13HRT4820	36
		0	3.5	Violet	70	AM1L18HRT4810	36
		1	4	Violet	150 loop	AM1L15HRT48L01	36
• RIBBED							
HRT 48 PRIME NEEDLE		2/0	3	Violet	90	AMP1L13HRT4820	36
		0	3.5	Violet	70	AMP1L18HRT4810	36
		1	4	Violet	150 loop	AMP1L15HRT48L01	36
• RIBBED							

Premisut PGLA - 910

▼ 1/2 circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HS 22		3/0	2	Violet	70	AM1L18HS2230	36
		3/0	2	Undyed	70	AM1L28HS2230	36
		2/0	3	Violet	70	AM1L18HS2220	36
		2/0	3	Undyed	75	AM1L22HS2220	36
• RIBBED • RTT							
HS 22 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18HS2230	36
		3/0	2	Undyed	70	AMP1L28HS2230	36
		2/0	3	Violet	70	AMP1L18HS2220	36
		2/0	3	Undyed	75	AMP1L22HS2220	36
• RIBBED • RTT							
HS 22 heavy PRIME NEEDLE		0	3.5	Undyed	70	AMP1L28HS2210H	36
		2	5	Violet	75	AMP1L12HS2202H	36
• RIBBED • RTT • HEAVY NEEDLE							
HS 26		3/0	2	Violet	70	AM1L18HS2630	36
		2/0	3	Undyed	70	AM1L28HS2620	36
		0	3.5	Violet	70	AM1L18HS2610	36
		1	4	Undyed	70	AM1L28HS2601	36
• RIBBED • RTT							
HS 26 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18HS2630	36
		2/0	3	Undyed	70	AMP1L28HS2620	36
		0	3.5	Violet	70	AMP1L18HS2610	36
		1	4	Undyed	70	AMP1L28HS2601	36
• RIBBED • RTT							

Premisut PGLA - 910

▼ 1/2 circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HS 36		3/0	2	Violet	70	AM1L18HS3630	36
		2/0	3	Violet	70	AM1L18HS3620	36
		2/0	3	Violet	90	AM1L13HS3620	36
		2/0	3	Undyed	90	AM1L23HS3620	36
		0	3.5	Violet	70	AM1L18HS3610	36
		0	3.5	Violet	90	AM1L13HS3610	36
		1	4	Violet	70	AM1L18HS3601	36
		1	4	Violet	90	AM1L13HS3601	36
		2	5	Violet	90	AM1L13HS3602*	36
• RIBBED • RTT *NO RTT							
HS 36 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18HS3630	36
		2/0	3	Violet	70	AMP1L18HS3620	36
		2/0	3	Violet	90	AMP1L13HS3620	36
		2/0	3	Undyed	90	AMP1L23HS3620	36
		0	3.5	Violet	70	AMP1L18HS3610	36
		0	3.5	Violet	90	AMP1L13HS3610	36
		1	4	Violet	70	AMP1L18HS3601	36
		1	4	Violet	90	AMP1L13HS3601	36
		2	5	Violet	90	AMP1L13HS3602*	36
• RIBBED • RTT *NO RTT							
HS 40		2/0	3	Violet	70	AM1L18HS4020	36
		0	3.5	Violet	70	AM1L18HS4010	36
		0	3.5	Violet	90	AM1L13HS4010	36
		1	4	Violet	70	AM1L18HS4001	36
		1	4	Violet	90	AM1L13HS4001	36
		2	5	Violet	90	AM1L13HS4002*	36
• RIBBED • RTT *NO RTT							
HS 40 PRIME NEEDLE		2/0	3	Violet	70	AMP1L18HS4020	36
		0	3.5	Violet	70	AMP1L18HS4010	36
		0	3.5	Violet	90	AMP1L13HS4010	36
		1	4	Violet	70	AMP1L18HS4001	36
		1	4	Violet	90	AMP1L13HS4001	36
		2	5	Violet	90	AMP1L13HS4002*	36
• RIBBED • RTT *NO RTT							

Premisut PGLA - 910

▼ 1/2 circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HS 48		0	3.5	Violet	90	AM1L13HS4810	36
		1	4	Violet	90	AM1L13HS4801	36
		2	5	Violet	90	AM1L13HS4802	36
• RIBBED							
HS 48 PRIME NEEDLE		0	3.5	Violet	90	AMP1L13HS4810	36
		1	4	Violet	90	AMP1L13HS4801	36
		2	5	Violet	90	AMP1L13HS4802	36
• RIBBED							

▼ Circle Extra Reverse Cutting – Long cutting edge (plastic/cosmetic)

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
XHS 16 PRIME NEEDLE		5/0	1	Undyed	45	AMP1L21XHS1650	36
		4/0	1.5	Undyed	45	AMP1L21XHS1640	36

▼ 3/8 Circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 12 PRIME NEEDLE		6/0	0.7	Violet	45	AMP1L11DS1260	12
		6/0	0.7	Undyed	45	AMP1L21DS1260	12
		5/0	1	Violet	45	AMP1L11DS1250	36
		4/0	1.5	Violet	45	AMP1L11DS1240	36

Premisut PGLA - 910

▼ 3/8 Circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 14 PRIME NEEDLE		6/0	0.7	Violet	45	AMP1L11DS1460	12
DS 16		6/0	0.7	Violet	45	AM1L11DS1660	12
		5/0	1	Violet	45	AM1L11DS1650	36
		5/0	1	Violet	75	AM1L12DS1650	36
		5/0	1	Undyed	75	AM1L22DS1650	36
		4/0	1.5	Violet	75	AM1L12DS1640	36
		4/0	1.5	Undyed	75	AM1L22DS1640	36
		3/0	2	Violet	75	AM1L12DS1630	36
DS 16 PRIME NEEDLE		6/0	0.7	Violet	45	AMP1L11DS1660	12
		5/0	1	Violet	45	AMP1L11DS1650	36
		5/0	1	Violet	75	AMP1L12DS1650	36
		5/0	1	Undyed	75	AMP1L22DS1650	36
		4/0	1.5	Violet	75	AMP1L12DS1640	36
		4/0	1.5	Undyed	75	AMP1L22DS1640	36
		3/0	2	Violet	75	AMP1L12DS1630	36
DS 19		5/0	1	Undyed	45	AM1L21DS1950	36
		5/0	1	Violet	45	AM1L11DS1950	36
		5/0	1	Violet	75	AM1L12DS1950	36
		4/0	1.5	Violet	45	AM1L11DS1940	36
		4/0	1.5	Violet	75	AM1L12DS1940	36
		4/0	1.5	Undyed	75	AM1L22DS1940	36
		3/0	2	Violet	75	AM1L12DS1930	36
		3/0	2	Undyed	75	AM1L22DS1930	36
• RIBBED • RTT							
DS 19 PRIME NEEDLE		5/0	1	Undyed	45	AMP1L21DS1950	36
		5/0	1	Violet	45	AMP1L11DS1950	36
		5/0	1	Violet	75	AMP1L12DS1950	36
		4/0	1.5	Violet	45	AMP1L11DS1940	36
		4/0	1.5	Violet	75	AMP1L12DS1940	36
		4/0	1.5	Undyed	75	AMP1L22DS1940	36
		3/0	2	Violet	75	AMP1L12DS1930	36
		3/0	2	Undyed	75	AMP1L22DS1930	36
• RIBBED • RTT							

Premisut PGLA - 910

3/8 Circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 19 fine needle PRIME NEEDLE		4/0	1.5	Violet	45	AMP1L11DS1940FN	36
		4/0	1.5	Undyed	45	AMP1L21DS1940FN	36
		4/0	1.5	Violet	70	AMP1L18DS1940FN	36
		4/0	1.5	Undyed	70	AMP1L28DS1940FN	36
• RTT							
DS 24		4/0	1.5	Violet	70	AM1L18DS2440	36
		4/0	1.5	Undyed	70	AM1L28DS2440	36
		3/0	2	Violet	75	AM1L12DS2430	36
		3/0	2	Undyed	75	AM1L22DS2430	36
		2/0	3	Violet	75	AM1L12DS2420	36
		2/0	3	Undyed	75	AM1L22DS2420	36
		0	3.5	Violet	75	AM1L12DS2410	36
• RIBBED • RTT							
DS 24 PRIME NEEDLE		4/0	1.5	Violet	70	AMP1L18DS2440	36
		4/0	1.5	Undyed	70	AMP1L28DS2440	36
		3/0	2	Violet	75	AMP1L12DS2430	36
		3/0	2	Undyed	75	AMP1L22DS2430	36
		2/0	3	Violet	75	AMP1L12DS2420	36
		2/0	3	Undyed	75	AMP1L22DS2420	36
		0	3.5	Violet	75	AMP1L12DS2410	36
• RIBBED • RTT							
DS 26		3/0	2	Violet	75	AM1L12DS2630	36
		2/0	3	Violet	75	AM1L12DS2620	36
• RIBBED • RTT							
DS 26 PRIME NEEDLE		3/0	2	Violet	75	AMP1L12DS2630	36
		2/0	3	Violet	75	AMP1L12DS2620	36
• RIBBED • RTT							

Premisut PGLA - 910

▼ 3/8 Circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 30		3/0	2	Violet	75	AM1L12DS3030	36
		2/0	3	Violet	75	AM1L12DS3020	36
		2/0	3	Undyed	75	AM1L22DS3020	36
		0	3.5	Violet	75	AM1L12DS3010	36
		0	3.5	Undyed	70	AM1L28DS3010	36
• RIBBED • RTT							
DS 30 PRIME NEEDLE		3/0	2	Violet	75	AMP1L12DS3030	36
		2/0	3	Violet	75	AMP1L12DS3020	36
		2/0	3	Undyed	75	AMP1L22DS3020	36
		0	3.5	Violet	75	AMP1L12DS3010	36
		0	3.5	Undyed	70	AMP1L28DS3010	36
• RIBBED • RTT							
DS 36		2/0	3	Violet	75	AM1L12DS3620	36
		0	3.5	Violet	75	AM1L12DS3610	36
		1	4	Violet	75	AM1L12DS3601	36
		2	5	Violet	75	AM1L12DS3602	36
• RIBBED • RTT							
DS 36 PRIME NEEDLE		2/0	3	Violet	75	AMP1L12DS3620	36
		0	3.5	Violet	75	AMP1L12DS3610	36
		1	4	Violet	75	AMP1L12DS3601	36
		2	5	Violet	75	AMP1L12DS3602	36
• RIBBED • RTT							
DS 40		0	3.5	Violet	75	AM1L12DS4010	36
		1	4	Violet	75	AM1L12DS4001	36
• RIBBED • RTT							
DS 40 PRIME NEEDLE		0	3.5	Violet	75	AMP1L12DS4010	36
		1	4	Violet	75	AMP1L12DS4001	36
• RIBBED • RTT							

Premisut PGLA - 910

▲ Circle Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
NS 16 PRIME NEEDLE		5/0	1	Violet	45	AMP1L11NS1650	36

▼ 3/8 Circle Reverse Extra Cutting – Long cutting edge

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2 x XDS 8 PRIME NEEDLE		6/0	0.7	Violet	45	AMP1L11XDS860-2	12
XDS 11 PRIME NEEDLE		6/0	0.7	Undyed	45	AMP1L21XDS1160	12
		5/0	1	Undyed	45	AMP1L21XDS1150	36
2 x XDS 11 PRIME NEEDLE		6/0	0.7	Violet	45	AMP1L11XDS1160-2	12
XDS 13 PRIME NEEDLE		6/0	0.7	Undyed	45	AMP1L21XDS1360	12
		5/0	1	Undyed	45	AMP1L21XDS1350	36
		4/0	1.5	Undyed	45	AMP1L21XDS1340	36

Premisut PGLA - 910

 **3/8 Circle Reverse Cutting – Long cutting edge (Extra Cutting)** 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2 x XDS 13 PRIME NEEDLE		5/0	1	Violet	45	AMP1L11XDS1350-2	36
XDS 16 PRIME NEEDLE		5/0	1	Undyed	45	AMP1L21XDS1650	36
		4/0	1.5	Undyed	45	AMP1L21XDS1640	36
XDS 19 PRIME NEEDLE		5/0	1	Undyed	45	AMP1L21XDS1950	36
		4/0	1.5	Undyed	45	AMP1L21XDS1940	36
		3/0	2	Undyed	45	AMP1L21XDS1930	36
• RTT							
XDS 24 PRIME NEEDLE		3/0	2	Undyed	45	AMP1L21XDS2430	36
		2/0	3	Undyed	45	AMP1L21XDS2420	36
• RTT							
XDS 26 PRIME NEEDLE		3/0	2	Undyed	70	AMP1L28XDS2630	36
• RTT							
XDS 30 PRIME NEEDLE		2/0	3	Violet	70	AMP1L18XDS3020	36
• RTT							

Premisut PGLA - 910

3/8 Circle Extra Cutting – Long cutting edge (plastic/cosmetic)

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
XNS 16 PRIME NEEDLE		5/0	1	Undyed	45	AMP1L21XNS1650	36
		4/0	1.5	Undyed	45	AMP1L21XNS1640	36
XNS 19 PRIME NEEDLE		3/0	2	Undyed	70	AMP1L28XNS1930	36
• RTT							
XNS 26 PRIME NEEDLE		4/0	1.5	Undyed	75	AMP1L22XNS2640	36
		3/0	2	Undyed	75	AMP1L22XNS2630	36
• RTT							

3/8 circle round bodied TAPERCUT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DRT 13 PRIME NEEDLE		6/0	0.7	Violet	70	AMP1L18DRT1360*	12
		5/0	1	Violet	70	AMP1L18DRT1350	36
DRT 17 PRIME NEEDLE		4/0	1.5	Undyed	70	AMP1L28DRT1740	36

Premisut PGLA - 910

● 3/8 Circle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DR 26 PRIME NEEDLE		3/0	2	Violet	75	AMP1L12DR2630	36
		2/0	3	Violet	75	AMP1L12DR2620	36
		• RIBBED • RTT					

○ 3/8 Circle Round bodied Blunt Point

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DRN 64		1	4	Violet	75	AM1L12DRN6401	36
		1	4	Violet	90	AM1L13DRN6401	36
		1	4	Violet	100	AM1L17DRN6401	36
		• RIBBED					
DRN 64 PRIME NEEDLE		1	4	Violet	75	AMP1L12DRN6401	36
		1	4	Violet	90	AMP1L13DRN6401	36
		1	4	Violet	100	AMP1L17DRN6401	36
		• RIBBED					

▲ Straight Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2xGS 13 PRIME NEEDLE		4/0	1.5	Violet	45	AMP1L11GS1340-2	36

Premisut PGLA - 910

▲ Straight Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
GS 25 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18GS2530	36
• RIBBED							
GS 48 PRIME NEEDLE		3/0	2	Undyed	70	AMP1L28GS4830	36
• RIBBED							
GS 60 PRIME NEEDLE		4/0	1.5	Undyed	70	AMP1L28GS6040	36
		3/0	2	Violet	75	AMP1L12GS6030	36
		3/0	2	Undyed	70	AMP1L28GS6030	36
		2/0	3	Violet	75	AMP1L12GS6020	36
		2/0	3	Undyed	70	AMP1L28GS6020	36
• RIBBED							

● Straight needle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
GR 19 PRIME NEEDLE		4/0	1.5	Violet	70	AMP1L18GR1940	36
		3/0	2	Violet	70	AMP1L18GR1930	36
• RIBBED							
2xGR 19 PRIME NEEDLE		4/0	1.5	Violet	45	AMP1L11GR1940-2	36
• RIBBED							

Premisut PGLA - 910

● Straight needle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
GR 48 PRIME NEEDLE		3/0	2	Violet	70	AMP1L18GR4830	36
• RIBBED							

● 5/8 Circle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
FR 27		3/0	2	Violet	75	AM1L12FR2730	36
		2/0	3	Violet	75	AM1L12FR2720	36
		0	3.5	Violet	70	AM1L18FR2710	36
		1	4	Violet	70	AM1L18FR2701	36
• RIBBED • RTT							
FR 27 PRIME NEEDLE		3/0	2	Violet	75	AMP1L12FR2730	36
		2/0	3	Violet	75	AMP1L12FR2720	36
		0	3.5	Violet	70	AMP1L18FR2710	36
		1	4	Violet	70	AMP1L18FR2701	36
• RIBBED • RTT							
2xFR 27		2/0	3	Violet	70	AM1L18FR2720-2	36
• RIBBED							
2xFR 27 PRIME NEEDLE		2/0	3	Violet	70	AMP1L18FR2720-2	36
• RIBBED							

Premisut PGLA - 910

● 5/8 Circle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
FR 36		2/0	3	Violet	70	AM1L18FR3620	36
		0	3.5	Violet	70	AM1L18FR3610	36
		1	4	Violet	70	AM1L18FR3601	36
FR 36 PRIME NEEDLE		2/0	3	Violet	70	AMP1L18FR3620	36
		0	3.5	Violet	70	AMP1L18FR3610	36
		1	4	Violet	70	AMP1L18FR3601	36
FR 40		2/0	3	Violet	70	AM1L18FR4020	36
		0	3.5	Violet	70	AM1L18FR4010	36
FR 40 PRIME NEEDLE		2/0	3	Violet	70	AMP1L18FR4020	36
		0	3.5	Violet	70	AMP1L18FR4010	36

● J-Hook Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
JR 30		2/0	3	Violet	70	AM1L18JR3020	36
		1	4	Violet	70	AM1L18JR3001	36

Premisut PGLA - 910

● J-Hook Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
JR 30 PRIME NEEDLE • RIBBED • RTT		2/0	3	Violet	70	AMP1L18JR3020	36
		1	4	Violet	70	AMP1L18JR3001	36
JR 36 • RIBBED • RTT		2/0	3	Violet	75	AM1L12JR3620	36
		0	3.5	Violet	75	AM1L12JR3610	36
JR 36 PRIME NEEDLE • RIBBED • RTT		2/0	3	Violet	75	AMP1L12JR3620	36
		0	3.5	Violet	75	AMP1L12JR3610	36

▼ J-Hook Round bodied TAPERCUT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
JRT 31 • RIBBED • RTT		1	4	Violet	75	AM1L12JRT3101	36
JRT 31 PRIME NEEDLE • RIBBED • RTT		1	4	Violet	75	AMP1L12JRT3101	36

Premisut PGLA - 910

1/2 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2xHSP 7 PRIME NEEDLE		8/0	0.4	Violet	30	AMP1L14HSP780-2	12
		7/0	0.5	Violet	45	AMP1L11HSP770-2	12
2xHSP 8 PRIME NEEDLE		6/0	0.7	Violet	45	AMP1L11HSP860-2	12
		6/0	0.7	Undyed	45	AMP1L21HSP860-2	12

3/8 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DSP 6.5 PRIME NEEDLE		8/0	0.4	Violet	30	AMP1L14DSP6.580	12
		7/0	0.5	Violet	30	AMP1L14DSP6.570	12
2xDSP 6.5 PRIME NEEDLE		8/0	0.4	Violet	30	AMP1L14DSP6.580-2	12
		7/0	0.5	Violet	45	AMP1L11DSP6.570-2	12
2xDSP 11 PRIME NEEDLE		5/0	1	Undyed	45	AMP1L21DSP1150-2	36

Premisut PGLA - 910

1/4 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2xVSP 6.5 PRIME NEEDLE		6/0	0.7	Violet	45	AMP1L11VSP6.560-2	12
2xVSP 8 PRIME NEEDLE		6/0	0.7	Violet	45	AMP1L11VSP860-2	12
		5/0	1	Violet	45	AMP1L11VSP850-2	36
		4/0	1.5	Violet	45	AMP1L11VSP840-2	36

Without needle

NEEDLE LENGTH MM	PRE-CUT SUTURES	USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
	without needle 180cm	4/0	1.5	Violet	180	AM1L1-18040	36
		3/0	2	Violet	180	AM1L1-18030	36
		2/0	3	Violet	180	AM1L1-18020	36
		0	3.5	Violet	180	AM1L1-18010	36
		1	4	Violet	180	AM1L1-18001	36
		2	5	Violet	180	AM1L1-18002	36
	without needle 150cm	4/0	1.5	Violet	150	AM1L1-15040	36
		3/0	2	Violet	150	AM1L1-15030	36
		2/0	3	Violet	150	AM1L1-15020	36
		0	3.5	Violet	150	AM1L1-15010	36
		1	4	Violet	150	AM1L1-15001	36
		2	5	Violet	150	AM1L1-15002	36
	without needle 3 x 50cm	3/0	2	Violet	3 x 50	AM1L1-3.5030	36
		2/0	3	Violet	3 x 50	AM1L1-3.5020	36
		0	3.5	Violet	3 x 50	AM1L1-3.5010	36
		1	4	Violet	3 x 50	AM1L1-3.5001	36
		2	5	Violet	3 x 50	AM1L1-3.5002	36

Premisut PGLA - 910

Without needle

NEEDLE LENGTH MM	PRE-CUT SUTURES	USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
	without needle 2 x 70cm	4/0	1.5	Violet	2 x 70	AM1L1-2.7040	36
		3/0	2	Violet	2 x 70	AM1L1-2.7030	36
		2/0	3	Violet	2 x 70	AM1L1-2.7020	36
		0	3.5	Violet	2 x 70	AM1L1-2.7010	36
		1	4	Violet	2 x 70	AM1L1-2.7001	36
		2	5	Violet	2 x 70	AM1L1-2.7002	36
	without needle 3 x 45cm	5/0	1	Violet	3 x 45	AM1L1-3.4550	36
		4/0	1.5	Violet	3 x 45	AM1L1-3.4540	36
		3/0	2	Violet	3 x 45	AM1L1-3.4530	36
		2/0	3	Violet	3 x 45	AM1L1-3.4520	36
		0	3.5	Violet	3 x 45	AM1L1-3.4510	36
		1	4	Violet	3 x 45	AM1L1-3.4501	36
	without needle 5 x 70cm	2	5	Violet	3 x 45	AM1L1-3.4502	36
		3/0	2	Violet	5 x 70	AM1L1-5.7030	36
		2/0	3	Violet	5 x 70	AM1L1-5.7020	36
		0	3.5	Violet	5 x 70	AM1L1-5.7010	36
		1	4	Violet	5 x 70	AM1L1-5.7001	36
		2	5	Violet	5 x 70	AM1L1-5.7002	36
	without needle 6 x 45cm	5/0	1	Violet	6 x 45	AM1L1-6.4550	36
		4/0	1.5	Violet	6 x 45	AM1L1-6.4540	36
		3/0	2	Violet	6 x 45	AM1L1-6.4530	36
		2/0	3	Violet	6 x 45	AM1L1-6.4520	36
		0	3.5	Violet	6 x 45	AM1L1-6.4510	36
		1	4	Violet	6 x 45	AM1L1-6.4501	36
	without needle 12 x 45cm	2	5	Violet	6 x 45	AM1L1-6.4502	36
		4/0	1.5	Violet	12 x 45	AM1L1-12.4540	36
		3/0	2	Violet	12 x 45	AM1L1-12.4530	36
		2/0	3	Violet	12 x 45	AM1L1-12.4520	36
		0	3.5	Violet	12 x 45	AM1L1-12.4510	36
		1	4	Violet	12 x 45	AM1L1-12.4501	36
		2	5	Violet	12 x 45	AM1L1-12.4502	36