

Description

Premisut SILK suture is a sterile multifilament braided non-absorbable surgical suture prepared from an organic protein called fibroin. This protein is derived from the domesticated species *Bombyx mori* of the family *Bombycidae*. Natural Silk contains approximately 30% extraneous material or gum which is removed by "Degumming" process. This process is essential for a compact braid while ensuring that the filaments retain their natural body and elasticity. Premisut SILK suture is dyed with black Hematoxylin and coated with bees wax.

Indications

Premisut SILK suture is intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological tissues.

Performance

Premisut SILK suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. While Premisut SILK suture is not absorbed, progressive degradation of the proteinaceous silk fibres in vivo may result in gradual loss of the suture's tensile strength over a period of time.

Tensile Strength

Exceeds USP & EP requirements

Total Absorbption

None

Color

Black



Sizes

8/0 to 2 USP

Usage Areas

General surgery, skin closure, ophthalmic, oral, gastrointestinal, gynecology, plastic, microsurgery, obstetrics.

Sterilization

Ethylene Oxide

Shelf Life

5 years

Packaging

Box of 12pcs / 36pcs



Premisut SILK

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 12 PRIME NEEDLE		6/0	0.7	Black	45	AMP541HR1260	12
HR 13		6/0	0.7	Black	45	AM541HR1360	12
		5/0	1	Black	75	AM542HR1350	36
		4/0	1.5	Black	45	AM541HR1340	36
HR 13 PRIME NEEDLE		6/0	0.7	Black	45	AMP541HR1360	12
		5/0	1	Black	75	AMP542HR1350	36
		4/0	1.5	Black	45	AMP541HR1340	36
HR 16		6/0	0.7	Black	45	AM541HR1660	12
		6/0	0.7	Black	75	AM542HR1660	12
		5/0	1	Black	75	AM542HR1650	36
		4/0	1.5	Black	75	AM542HR1640	36
		3/0	2	Black	75	AM542HR1630	36
		2/0	3	Black	75	AM542HR1620	36
HR 16 PRIME NEEDLE		6/0	0.7	Black	45	AMP541HR1660	12
		5/0	1	Black	75	AMP542HR1650	36
		4/0	1.5	Black	75	AMP542HR1640	36
		3/0	2	Black	75	AMP542HR1630	36
		2/0	3	Black	75	AMP542HR1620	36

Premisut SILK

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM			USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 20			5/0	1	Black	75	AM542HR2050	36
			4/0	1.5	Black	75	AM542HR2040	36
			3/0	2	Black	75	AM542HR2030	36
			2/0	3	Black	75	AM542HR2020	36
			0	3.5	Black	75	AM542HR2010	36
• RIBBED • RTT								
HR 20 PRIME NEEDLE			5/0	1	Black	75	AMP542HR2050	36
			4/0	1.5	Black	75	AMP542HR2040	36
			3/0	2	Black	75	AMP542HR2030	36
			2/0	3	Black	75	AMP542HR2020	36
			0	3.5	Black	75	AMP542HR2010	36
• RIBBED • RTT								
HR 24			3/0	2	Black	75	AM542HR2430	36
			2/0	3	Black	75	AM542HR2420	36
• RIBBED • RTT								
HR 24 PRIME NEEDLE			3/0	2	Black	75	AMP542HR2430	36
			2/0	3	Black	75	AMP542HR2420	36
• RIBBED • RTT								
HR 26			4/0	1.5	Black	75	AM542HR2640	36
			3/0	2	Black	75	AM542HR2630	36
			2/0	3	Black	75	AM542HR2620	36
			0	3.5	Black	75	AM542HR2610	36
			1	4	Black	75	AM542HR2601	36
• RIBBED • RTT								

Premisut SILK

● ½ Circle Round bodied TAPERPOINT 

NEEDLE LENGTH MM			USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 26 PRIME NEEDLE			4/0	1.5	Black	75	AMP542HR2640	36
			3/0	2	Black	75	AMP542HR2630	36
			2/0	3	Black	75	AMP542HR2620	36
			0	3.5	Black	75	AMP542HR2610	36
			1	4	Black	75	AMP542HR2601	36
• RIBBED • RTT								
HR 30			3/0	2	Black	75	AM542HR3030	36
			2/0	3	Black	75	AM542HR3020	36
			0	3.5	Black	75	AM542HR3010	36
			1	4	Black	75	AM542HR3001	36
• RIBBED • RTT								
HR 30 PRIME NEEDLE			3/0	2	Black	75	AMP542HR3030	36
			2/0	3	Black	75	AMP542HR3020	36
			0	3.5	Black	75	AMP542HR3010	36
			1	4	Black	75	AMP542HR3001	36
• RIBBED • RTT								
HR 37			2/0	3	Black	75	AM542HR3720	36
			0	3.5	Black	75	AM542HR3710	36
			1	4	Black	75	AM542HR3701	36
			2	5	Black	75	AM542HR3702	36
• RIBBED • RTT								
HR 37 PRIME NEEDLE			2/0	3	Black	75	AMP542HR3720	36
			0	3.5	Black	75	AMP542HR3710	36
			1	4	Black	75	AMP542HR3701	36
			2	5	Black	75	AMP542HR3702	36
• RIBBED • RTT								

Premisut SILK

● ½ Circle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HR 40		2/0	3	Black	75	AM542HR4020	36
		0	3.5	Black	75	AM542HR4010	36
		1	4	Black	75	AM542HR4001	36
		2	5	Black	75	AM542HR4002	36
• RIBBED • RTT							
HR 40 PRIME NEEDLE		2/0	3	Black	75	AMP542HR4020	36
		0	3.5	Black	75	AMP542HR4010	36
		1	4	Black	75	AMP542HR4001	36
		2	5	Black	75	AMP542HR4002	36
• RIBBED • RTT							

▼ 1/2 circle round bodied TAPERCUT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRT 18		4/0	1.5	Black	75	AM542HRT1840	36
		3/0	2	Black	45	AM541HRT1830	36
		3/0	2	Black	75	AM542HRT1830	36
HRT 18 PRIME NEEDLE		4/0	1.5	Black	75	AMP542HRT1840	36
		3/0	2	Black	45	AMP541HRT1830	36
		3/0	2	Black	75	AMP542HRT1830	36

Premisut SILK

 **1/2 circle round bodied TAPERCUT** 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HRT 20 PRIME NEEDLE		4/0	1.5	Black	75	AMP542HRT2040	36
• RIBBED • RTT							
HRT 26		3/0	2	Black	75	AM542HRT2630	36
• RIBBED • RTT		2/0	3	Black	75	AM542HRT2620	36
HRT 26 PRIME NEEDLE		3/0	2	Black	75	AMP542HRT2630	36
• RIBBED • RTT		2/0	3	Black	75	AMP542HRT2620	36
HRT 30		3/0	2	Black	75	AM542HRT3030	36
• RIBBED • RTT							

Premisut SILK

▼ 1/2 circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HS 13		5/0	1	Black	45	AM541HS1350	36
HS 13 PRIME NEEDLE		5/0	1	Black	45	AMP541HS1350	36
HS 17		4/0	1.5	Black	45	AM541HS1740	36
		4/0	1.5	Black	75	AM542HS1740	36
		3/0	2	Black	75	AM542HS1730	36
HS 17 PRIME NEEDLE		4/0	1.5	Black	45	AMP541HS1740	36
		4/0	1.5	Black	75	AMP542HS1740	36
		3/0	2	Black	75	AMP542HS1730	36
HS 22		3/0	2	Black	45	AM541HS2230	36
		3/0	2	Black	75	AM542HS2230	36
		2/0	3	Black	45	AM541HS2220	36
		2/0	3	Black	75	AM542HS2220	36
• RIBBED • RTT							

Premisut SILK

▼ 1/2 circle Reverse Cutting 

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HS 22 PRIME NEEDLE		3/0	2	Black	45	AMP541HS2230	36
		3/0	2	Black	75	AMP542HS2230	36
		2/0	3	Black	45	AMP541HS2220	36
		2/0	3	Black	75	AMP542HS2220	36
HS 26		3/0	2	Black	75	AM542HS2630	36
		2/0	3	Black	75	AM542HS2620	36
		0	3.5	Black	75	AM542HS2610	36
HS 26 PRIME NEEDLE		3/0	2	Black	75	AMP542HS2630	36
		2/0	3	Black	75	AMP542HS2620	36
HS 36		2/0	3	Black	75	AM542HS3620	36
		0	3.5	Black	75	AM542HS3610	36
		1	4	Black	75	AM542HS3601	36
HS 36 PRIME NEEDLE		2/0	3	Black	75	AMP542HS3620	36
		0	3.5	Black	75	AMP542HS3610	36
		1	4	Black	75	AMP542HS3601	36

Premisut SILK

▼ 1/2 circle Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
HNS 26		3/0	2	Black	75	AM542HNS2630	36
		2/0	3	Black	75	AM542HNS2620	36
• RTT							
HNS 30		2/0	3	Black	75	AM542HNS3020	36
• RTT							

● 3/8 Circle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DR 9.3 PRIME NEEDLE		7/0	0.5	Black	30	AMP544DR9.370	12
DR 13 PRIME NEEDLE		6/0	0.7	Black	75	AMP542DR1360	12
		5/0	1	Black	75	AMP542DR1350	36
		4/0	1.5	Black	75	AMP542DR1340	36
2 x DR 13 PRIME NEEDLE		6/0	0.5	Black	45	AMP541DR1360-2	12

Premisut SILK

● 3/8 Circle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DR 17 PRIME NEEDLE		4/0	1	Black	75	AMP542DR1740	36

▼ 3/8 circle round bodied TAPERCUT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DRT 18 PRIME NEEDLE		4/0	1.5	Black	75	AMP542DRT1840	36
		3/0	2	Black	75	AMP542DRT1830	36

▼ 3/8 Circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 8 PRIME NEEDLE		7/0	0.5	Black	30	AMP544DS870	12
		6/0	0.7	Black	30	AMP544DS860	12
2 x DS 8 PRIME NEEDLE		7/0	0.5	Black	45	AMP541DS870-2	12
		6/0	0.7	Black	45	AMP541DS860-2	12

Premisut SILK

▼ 3/8 Circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 12 PRIME NEEDLE		7/0	0.5	Black	45	AMP541DS1270	12
		6/0	0.7	Black	45	AMP541DS1260	12
		5/0	1	Black	45	AMP541DS1250	36
		4/0	1.5	Black	45	AMP541DS1240	36
		3/0	2	Black	75	AMP542DS1230	36
DS 14		6/0	0.7	Black	45	AM541DS1460	12
		5/0	1	Black	45	AM541DS1450	36
		4/0	1.5	Black	45	AM541DS1440	36
DS 14 PRIME NEEDLE		6/0	0.7	Black	45	AMP541DS1460	12
		5/0	1	Black	45	AMP541DS1450	36
		4/0	1.5	Black	45	AMP541DS1440	36
DS 16		6/0	0.7	Black	45	AM541DS1660	12
		6/0	0.7	Black	75	AM542DS1660	12
		5/0	1	Black	75	AM542DS1650	36
		4/0	1.5	Black	75	AM542DS1640	36
		3/0	2	Black	75	AM542DS1630	36
		2/0	3	Black	75	AM542DS1620	36
DS 16 PRIME NEEDLE		6/0	0.7	Black	45	AMP541DS1660	12
		6/0	0.7	Black	75	AMP542DS1660	12
		5/0	1	Black	75	AMP542DS1650	36
		4/0	1.5	Black	75	AMP542DS1640	36
		3/0	2	Black	75	AMP542DS1630	36
		2/0	3	Black	75	AMP542DS1620	36

Premisut SILK

▼ 3/8 Circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 19		6/0	0.7	Black	45	AM541DS1960*	12
		6/0	0.7	Black	75	AM542DS1960*	12
		5/0	1	Black	45	AM541DS1950	36
		5/0	1	Black	75	AM542DS1950	36
		4/0	1.5	Black	45	AM541DS1940	36
		4/0	1.5	Black	75	AM542DS1940	36
		3/0	2	Black	45	AM541DS1930	36
		3/0	2	Black	75	AM542DS1930	36
		2/0	3	Black	45	AM541DS1920	36
		2/0	3	Black	75	AM542DS1920	36
		0	3.5	Black	75	AM542DS1910	36
		1	4	Black	75	AM542DS1901	36
		• RIBBED • RTT * NO RTT					
DS 19 PRIME NEEDLE		6/0	0.7	Black	45	AMP541DS1960*	12
		6/0	0.7	Black	75	AMP542DS1960*	12
		5/0	1	Black	45	AMP541DS1950	36
		5/0	1	Black	75	AMP542DS1950	36
		4/0	1.5	Black	45	AMP541DS1940	36
		4/0	1.5	Black	75	AMP542DS1940	36
		3/0	2	Black	45	AMP541DS1930	36
		3/0	2	Black	75	AMP542DS1930	36
		2/0	3	Black	45	AMP541DS1920	36
		2/0	3	Black	75	AMP542DS1920	36
		0	3.5	Black	75	AMP542DS1910	36
		1	4	Black	75	AMP542DS1901	36
		• RIBBED • RTT * NO RTT					
DS 24		6/0	0.7	Black	75	AM542DS2460*	12
		5/0	1	Black	75	AM542DS2450	36
		4/0	1.5	Black	45	AM541DS2440	36
		4/0	1.5	Black	75	AM542DS2440	36
		3/0	2	Black	45	AM541DS2430	36
		3/0	2	Black	75	AM542DS2430	36
		2/0	3	Black	75	AM542DS2420	36
		0	3.5	Black	75	AM542DS2410	36
		1	4	Black	75	AM542DS2401	36
		• RIBBED • RTT *NO RTT					

Premisut SILK

▼ 3/8 Circle Reverse Cutting

NEEDLE LENGTH MM			USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 24 PRIME NEEDLE			6/0	0.7	Black	75	AMP542DS2460*	12
			5/0	1	Black	75	AMP542DS2450	36
			4/0	1.5	Black	45	AMP541DS2440	36
			4/0	1.5	Black	75	AMP542DS2440	36
			3/0	2	Black	45	AMP541DS2430	36
			3/0	2	Black	75	AMP542DS2430	36
			2/0	3	Black	75	AMP542DS2420	36
			0	3.5	Black	75	AMP542DS2410	36
			1	4	Black	75	AMP542DS2401	36
• RIBBED • RTT *NO RTT								
DS 26			3/0	2	Black	75	AM542DS2630	36
			2/0	3	Black	75	AM542DS2620	36
			0	3.5	Black	75	AM542DS2610	36
			1	4	Black	75	AM542DS2601	36
• RIBBED • RTT								
DS 26 PRIME NEEDLE			3/0	2	Black	75	AMP542DS2630	36
			2/0	3	Black	75	AMP542DS2620	36
			0	3.5	Black	75	AMP542DS2610	36
			1	4	Black	75	AMP542DS2601	36
• RIBBED • RTT								
DS 30			5/0	1	Black	75	AM542DS3050	36
			4/0	1.5	Black	75	AM542DS3040	36
			3/0	2	Black	75	AM542DS3030	36
			2/0	3	Black	75	AM542DS3020	36
			0	3.5	Black	75	AM542DS3010	36
			1	4	Black	75	AM542DS3001	36
• RIBBED • RTT								

Premisut SILK

▼ 3/8 Circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 30 PRIME NEEDLE		5/0	1	Black	75	AMP542DS3050	36
		4/0	1.5	Black	75	AMP542DS3040	36
		3/0	2	Black	75	AMP542DS3030	36
		2/0	3	Black	75	AMP542DS3020	36
		0	3.5	Black	75	AMP542DS3010	36
		1	4	Black	75	AMP542DS3001	36
DS 37		3/0	2	Black	75	AM542DS3730	36
		2/0	3	Black	75	AM542DS3720	36
		0	3.5	Black	75	AM542DS3710	36
		1	4	Black	75	AM542DS3701	36
		2	5	Black	75	AM542DS3702	36
DS 37 PRIME NEEDLE		3/0	2	Black	75	AMP542DS3730	36
		2/0	3	Black	75	AMP542DS3720	36
		0	3.5	Black	75	AMP542DS3710	36
		1	4	Black	75	AMP542DS3701	36
		2	5	Black	75	AMP542DS3702	36
DS 40		2/0	3	Black	75	AM542DS4020	36
		0	3.5	Black	75	AM542DS4010	36
		1	2	Black	75	AM542DS4001	36
		2	5	Black	90	AM543DS4002*	36
DS 40 PRIME NEEDLE		2/0	3	Black	75	AMP542DS4020	36
		0	3.5	Black	75	AMP542DS4010	36
		1	2	Black	75	AMP542DS4001	36
		2	5	Black	90	AMP543DS4002*	36

Premisut SILK

▼ 3/8 Circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 45		2/0	3	Black	75	AM542DS4520	36
		0	3.5	Black	75	AM542DS4510	36
		1	4	Black	75	AM542DS4501	36
		2	5	Black	75	AM542DS4502	36
• RIBBED • RTT							
DS 45 PRIME NEEDLE		2/0	3	Black	75	AMP542DS4520	36
		0	3.5	Black	75	AMP542DS4510	36
		1	4	Black	75	AMP542DS4501	36
		2	5	Black	75	AMP542DS4502	36
• RIBBED • RTT							
DS 50		2/0	3	Black	90	AM543DS5020	36
		0	3.5	Black	90	AM543DS5010	36
		1	4	Black	90	AM543DS5001	36
• RIBBED							
DS 50 PRIME NEEDLE		2/0	3	Black	90	AMP543DS5020	36
		0	3.5	Black	90	AMP543DS5010	36
		1	4	Black	90	AMP543DS5001	36
• RIBBED							
DS 60		2/0	3	Black	75	AM542DS6020	36
		0	3.5	Black	75	AM542DS6010	36
		1	4	Black	75	AM542DS6001	36
• RIBBED							
DS 60 PRIME NEEDLE		2/0	3	Black	75	AMP542DS6020	36
		0	3.5	Black	75	AMP542DS6010	36
		1	4	Black	75	AMP542DS6001	36
• RIBBED							

Premisut SILK

▼ 3/8 Circle Reverse Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
DS 90		1	4	Black	90	AM543DS9001	36
		2	5	Black	90	AM543DS9002	36
• RIBBED							
DS 90 PRIME NEEDLE		1	4	Black	90	AMP543DS9001	36
		2	5	Black	90	AMP543DS9002	36
• RIBBED							

▲ 3/8 Circle Conventional Extra Cutting – Long cutting edge (Extra Cutting/cosmetic)

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
XNS 16 PRIME NEEDLE		6/0	0.7	Black	45	AMP541XNS1660	12
		5/0	1	Black	45	AMP541XNS1650	36
		4/0	1.5	Black	45	AMP541XNS1640	36

▼ 3/8 Circle Reverse Cutting – Long cutting edge (Extra Cutting/cosmetic)

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
XDS 11 PRIME NEEDLE		6/0	0.7	Black	45	AMP541XDS1160	12

Premisut SILK

3/8 Circle Reverse Cutting – Long cutting edge (Extra Cutting/cosmetic)

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2 x XDS 11 PRIME NEEDLE		6/0	0.7	Black	45	AMP541XDS1160-2	12
XDS 13 PRIME NEEDLE		5/0	1	Black	45	AMP541XDS1350	36
		4/0	1.5	Black	45	AMP541XDS1340	36

● Straight needle Round bodied TAPERPOINT

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
GR 19 PRIME NEEDLE		4/0	1.5	Black	45	AMP541GR1940	36
		4/0	1.5	Black	75	AMP542GR1940	36
• RIBBED							
GR 40 PRIME NEEDLE		2/0	3	Black	75	AMP542GR4020	36
• RIBBED							

Premisut SILK

▲ Straight Cutting

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
GS 40		2/0	3	Black	75	AM542GS4020	36
• RIBBED							
GS 40 PRIME NEEDLE		2/0	3	Black	75	AMP542GS4020	36
• RIBBED							
GS 55		2/0	3	Black	75	AM542GS5520	36
• RIBBED							
GS 60		3/0	2	Black	75	AM542GS6030	36
		2/0	3	Black	75	AM542GS6020	36
		0	3.5	Black	75	AM542GS6010	36
		1	4	Black	75	AM542GS6001	36
• RIBBED							
GS 60 PRIME NEEDLE		3/0	2	Black	75	AMP542GS6030	36
		2/0	3	Black	75	AMP542GS6020	36
		0	3.5	Black	75	AMP542GS6010	36
		1	4	Black	75	AMP542GS6001	36
• RIBBED							

Premisut SILK

3/8 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2 x DSP 6.5 PRIME NEEDLE		8/0	0.4	Black	45	AMP541DSP6.580-2	12
		7/0	0.5	Black	45	AMP541DSP6.570-2	12
		6/0	0.7	Black	45	AMP541DSP6.560-2	12

1/4 Circle MICROPOINT SPATULA

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
2 x VSP 8 PRIME NEEDLE		6/0	0.7	Black	45	AMP541VSP860-2	12
		5/0	1	Black	45	AMP541VSP850-2	36

Without needle

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
	without needle 150cm	3/0	2	Black	150	AM54-15030	36
		2/0	3	Black	150	AM54-15020	36
		0	3.5	Black	150	AM54-15010	36
		1	4	Black	150	AM54-15001	36
		2	5	Black	150	AM54-15002	36
		3&4	6	Black	150	AM54-15004	36
	without needle 180cm	3/0	2	Black	180	AM54-18030	36
		2/0	3	Black	180	AM54-18020	36
		0	3.5	Black	180	AM54-18010	36
		1	4	Black	180	AM54-18001	36
		2	5	Black	180	AM54-18002	36
		3&4	6	Black	180	AM54-18004	36
	without needle 2 x 70cm	3/0	2	Black	2 x 70	AM54-2.7030	36
		2/0	3	Black	2 x 70	AM54-2.7020	36
		0	3.5	Black	2 x 70	AM54-2.7010	36
		1	4	Black	2 x 70	AM54-2.7001	36
		2	5	Black	2 x 70	AM54-2.7002	36

Premisut SILK

Without needle

NEEDLE LENGTH MM		USP	EP	COLOUR	THREAD LENGTH CM	ATLAS REF	PCS/BOX
	without needle 6 x 45cm	4/0	1.5	Black	6 x 45	AM54-6.4540	36
		3/0	2	Black	6 x 45	AM54-6.4530	36
		2/0	3	Black	6 x 45	AM54-6.4520	36
		0	3.5	Black	6 x 45	AM54-6.4510	36
		1	4	Black	6 x 45	AM54-6.4501	36
		2	5	Black	6 x 45	AM54-6.4502	36
	without needle 17 x 45cm	4/0	1.5	Black	17 x 45	AM54-17.4540	36
		3/0	2	Black	17 x 45	AM54-17.4530	36
		2/0	3	Black	17 x 45	AM54-17.4520	36
		0	3.5	Black	17 x 45	AM54-17.4510	36
		1	4	Black	17 x 45	AM54-17.4501	36
		2	5	Black	17 x 45	AM54-17.4502	36
	without needle 13 x 60cm	4/0	1.5	Black	13 x 60	AM54-13.6040	36
		3/0	2	Black	13 x 60	AM54-13.6030	36
		2/0	3	Black	13 x 60	AM54-13.6020	36
		0	3.5	Black	13 x 60	AM54-13.6010	36
		1	4	Black	13 x 60	AM54-13.6001	36
		2	5	Black	13 x 60	AM54-13.6002	36
	without needle 10 x 75cm	4/0	1.5	Black	10 x 75	AM54-10.7540	36
		3/0	2	Black	10 x 75	AM54-10.7530	36
		2/0	3	Black	10 x 75	AM54-10.7520	36
		0	3.5	Black	10 x 75	AM54-10.7510	36
		1	4	Black	10 x 75	AM54-10.7501	36
		2	5	Black	10 x 75	AM54-10.7502	36