

Description

Premisut UCT cotton tape is a white woven cotton ligature.

Indications

Premisut UCT cotton tape is indicated for use in tying the umbilical cord of new born infant. It is also frequently used in pediatric and cardiovascular operative procedures to suspend small structures and vessels. Umbilical tape easily absorbs blood when used in an area of gross bleeding.

Performance

Premisut UCT cotton tape is strong enough to tie off the umbilical cord and for suspension of small structures and vessels. It elicits a minimal initial inflammatory reaction in tissues.

Tensile Strength

Permanent (Exceeds USP & EP requirements)

Total Absorption

Permanent

Color

White Woven cotton ligature

Usage Areas

Tying the umbilical cord of new born infant, pediatric and cardiovascular operative procedures to suspend small structures and vessels.

Sterilization

Gamma

Shelf Life

5 years

Packaging

Box of 12pcs

Ref	WIDTH	LENGTH	PCS/BOX
AMUCT45	3mm	45cm	12
AMUCT75	3mm	75CM	12